

amoena

Breast Care. Empowering women.



**PRODUCT CATALOGUE
2024/25**

Amoena. Your solution.

AMOENA BALANCE

PORTFOLIO OVERVIEW

Amoena breast shapers balance out asymmetry after breast-conserving therapy, breast reconstruction, multiple biopsies, or in the case of natural uneven breast development. The soft silicone shapers fit snugly against the natural curves of the remaining breast tissue and fill the bra cup for a symmetrical silhouette, preventing tension, improving posture and correcting muscular imbalances.

AVAILABLE COLOURS



- **Invisible** – fine edges ensure a natural fit between the shaper and the body
- **Versatile** – various shapes, sizes and thicknesses available for optimal customisation
- **Simple** – required thickness can be determined using our simple, tried-and-tested cup difference method

BALANCE ADAPT AIR

Breast shaper with an adjustable back (page 34)



BALANCE CONTACT

Breast shaper with a self-adhesive back (page 34)



BALANCE NATURA

Breast shaper with a temperature-regulating back (page 35)



BALANCE ESSENTIAL

Breast shaper for basic care (page 35)



BALANCE – BREAST SHAPERS



BALANCE ADAPT AIR

Breast shaper with an adjustable back

- Precisely and individually adaptable thanks to integrated AirChamber Technology – the perfect fit every time, even for patients who experience weight fluctuations or swelling due to lymphoedema.
- Sits close to the body thanks to flexible, adjustable back
- Fine volume adjustment offers greater flexibility when choosing a bra
- Reduces heat and perspiration thanks to Comfort+ technology



Image	Product name	Item no.	Shape	Thickness	Suitable for cup difference	Sizes	Colour	Layers
	Balance Adapt Air MD	233	Delta	Average/volume (adjustable)	1–2cm and >2 cm	2–12	Ivory	3
	Balance Adapt Air SV	234	Varia	Special shape	Cup difference method not applicable	S–XL	Ivory	3

BALANCE CONTACT

Breast shaper with a self-adhesive back

- Adapts especially well to the movement of the body
- Stays firmly in place
- Reduces heat and perspiration thanks to Comfort+ technology



Image	Product name	Item no.	Shape	Thickness	Suitable for cup difference	Sizes	Colour	Layers
	Balance Contact MD	229	Delta	Medium	1–2cm	2–12	Ivory	2
	Balance Contact VD	230	Delta	Volume	> 2cm	3/4–9/10	Ivory	2
	Balance Contact SV	286	Varia	Special shape	Cup difference method not applicable	S–XL	Ivory	2

BALANCE NATURA

Breast shaper with a temperature-regulating back

- Reduces heat and perspiration thanks to Comfort+ technology

Advanced



Image	Product name	Item no.	Shape	Thickness	Suitable for cup difference	Sizes	Colour	Layers
	Balance Natura TO	227	Oval	Thin	< 1cm	S–XL	Ivory	2
	Balance Natura TD	217	Delta	Thin	< 1cm	1–10	Ivory	2
	Balance Natura MD	220	Delta	Medium	1–2cm	2–12	Ivory	2
	Balance Natura VD	222	Delta	Volume	> 2cm	3/4–9/10	Ivory	2
	Balance Natura Light VD	221	Delta	Volume, 25% lighter	> 2cm	3/4–11/12	Ivory	2
	Balance Natura SE	231	Ellipse	Special shape	Cup difference method not applicable	S–XL	Ivory	2
	Balance Natura SV	285	Varia	Special shape	Cup difference method not applicable	S–XL	Ivory	2

Basic



BALANCE ESSENTIAL

Breast shaper that provides the basic care needed for a symmetrical appearance

Image	Product name	Item no.	Shape	Thickness	Suitable for cup difference	Sizes	Colour	Layers
	Balance Essential TO	228	Oval	Thin	< 1cm	S–XL	Ivory	1
	Balance Essential TD	218	Delta	Thin	< 1cm	1–10	Ivory	1
	Balance Essential MD	223	Delta	Average	1–2cm	2–12	Ivory Tawny 4–10	1
	Balance Essential VD	225	Delta	Volume	> 2cm	3/4–9/10	Ivory	1
	Balance Essential Light VD	224	Delta	Volume, 25% lighter	> 2cm	3/4–11/12	Ivory Tawny 3/4 – 9/10	1
	Balance Essential SE	232	Ellipse	Special shape	Cup difference method not applicable	S–XL	Ivory	1

Sizing chart, fitting tips and advice regarding the cup difference method on page 92 and 93.

BREAST FORMS

PORTFOLIO OVERVIEW

AVAILABLE COLOURS



ADAPT AIR

Breast form with an adjustable back
(page 40)



Premium



CONTACT

Breast form with a self-adhesive back
(page 42)




Premium



ENERGY

Breast form with a breathable back
(page 44)



Advanced



NATURA

Breast form with a temperature-regulating back
(page 46)



Advanced



ESSENTIAL

Breast form that provides the basic care needed for a symmetrical appearance
(page 48)



Basic



ADAPT AIR

Breast form with an adjustable back



ADAPT AIR triple-layer breast form

-  Soft, matte PU film and premium-quality silicone guarantee a pleasant feel
-  **Comfort+**
Reduces heat and perspiration behind the breast form
-  **AirChamber Technology**
 - Precisely and individually adaptable thanks to integrated AirChamber Technology – the perfect fit every time, even for patients who experience weight fluctuations or swelling due to lymphoedema
 - Fine volume adjustment offers greater flexibility when choosing a bra
 - Reduces heat and perspiration thanks to Comfort+ material
 - Lightweight: InTouch Light and Xtra Light silicone make the breast form 25-40% lighter
-  **Flexible, individually adjustable back**
Particularly close fit even with uneven chest wall.

PARTICULARLY SUITABLE FOR WOMEN WHO:

- have uneven scar tissue
- experience weight fluctuations (and the resulting fluctuations in cup size)
- have swelling from lymphoedema



WEIGHT OPTIONS



COLOURS



Premium

Overview

Image	Product name	Item no.	Cup form ¹	Weight	Shape ²	Sizes	Layers	Colour*
	Adapt Air Xtra Light 1S^N	328	1 – shallow	– 40%	S	1–14	3	Ivory
	Adapt Air Xtra Light 2S^N	326	2 – average	– 40%	S	1–14	3	Ivory Tawny (upon request)
	Adapt Air Light 1S^N	329	1 – shallow	– 25%	S	1–14	3	Ivory
	Adapt Air Light 2S^N	327	2 – average	– 25%	S	1–14	3	Ivory

For fitting tips and sizing chart see page 94 and 95.

1) OVERVIEW OF CUP FORMS:

1 - Shallow cup form: concave silhouette, natural shape with additional volume in the lower cup area and a flatter upper breast; wider basic form

2 – Average cup form: even silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a larger underbust and small to medium cup size

3 – Full cup form: convex silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a smaller underbust and large cup size

2) OVERVIEW OF FORMS:

S = Symmetrical: Can be worn on the left and right side interchangeably


A = Asymmetrical: Specifically designed to be worn on the left or right side respectively. Extended form to compensate for missing tissue under the arm

E = Extra: Specifically designed to be worn on the left or right side respectively. Particularly suitable for women with a smaller cup size and fuller breast shape

*Light und Xtra Light forms have a more muted colour.

CONTACT

Breast form with a self-adhesive back



CONTACT triple-layer breast form



Soft, matte PU film and premium-quality silicone guarantee a pleasant feel



Comfort+

Reduces heat and perspiration behind the breast form



Contact

- Adapts to body movements just like a natural breast
- Provides relief for neck and shoulders thanks to direct adhesion to the chest wall
- Offers more flexibility when choosing clothes

PARTICULARLY SUITABLE FOR WOMEN WHO:

- suffer from shoulder/neck pain
- have a larger cup size
- have lymphoedema

WEIGHT OPTIONS



Standard weight



-25%

COLOURS



Ivory



Tawny

Image	Product name	Item no.	Cup form ¹	Weight	Shape ²	Sizes	Colour*	Layers
	Contact Light 2S	380	2 – average	– 25%	S	1–14	Ivory	3
	Contact Light 3S	385	3 – volume	– 25%	S	0–10	Ivory	3
	Contact 1S	384	1 – shallow	full even weight distribution	S	1–14	Ivory	3
	Contact 2S	381	2 – average	full even weight distribution	S	1–12	Ivory	3
	Contact 2A	383	2 – average	full even weight distribution	A	0–12 L/R	Ivory	3
	Contact 3S	382	3 – volume	full even weight distribution	S	0–10	Ivory Tawny	3
	Contact 3E	386	3 – volume	full even weight distribution	E	1–14 L/R	Ivory Tawny (upon request)	3

For fitting tips and sizing chart see page 94 and 95.

1) OVERVIEW OF CUP FORMS:

1 - Shallow cup form: concave silhouette, natural shape with additional volume in the lower cup area and a flatter upper breast; wider basic form

2 – Average cup form: even silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a larger underbust and small to medium cup size

3 – Full cup form: convex silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a smaller underbust and large cup size

2) OVERVIEW OF FORMS:

S = Symmetrical: Can be worn on the left and right side interchangeably



A = Asymmetrical: Specifically designed to be worn on the left or right side respectively. Extended form to compensate for missing tissue under the arm

E = Extra: Specifically designed to be worn on the left or right side respectively. Particularly suitable for women with a smaller cup size and fuller breast shape

*Light und Xtra Light forms have a more muted colour.

ENERGY

Breast form with a breathable back



30% more Comfort+ than other breast forms for even better temperature regulation

Energy triple-layer breast form

- Soft, matte PU film and premium-quality silicone guarantee a pleasant feel
- Comfort+**
Reduces heat and perspiration behind the breast form
- 3D pearls**
Ensure optimal air circulation and allow moisture to evaporate

Energy Cosmetic triple-layer breast form

- Comfort+**
Reduces heat and perspiration behind the breast form
- +**
- Cosmetic**
Lightweight form with a natural appearance

PARTICULARLY SUITABLE FOR WOMEN WHO:

- experience hot flushes due to menopause or medication
- have a tendency to perspire in hotter temperatures
- take part in sports

WEIGHT OPTIONS



COLOURS



Advanced

Overview

Image	Product name	Item no.	Cup form ¹	Weight	Shape ²	Sizes	Colour*	Layers
	Energy Light 2 S	342	2 – average	– 25%	S	1–14	Ivory Tawny (upon request)	3
	Energy Cosmetic 2S*	310	2 – average	– 20%	S	1–14	Ivory	3
	Energy Cosmetic 3S*	311	3 – volume	– 20%	S	1–12	Ivory	3
	Energy 1S	349	1 – shallow	full even weight distribution	S	1–12	Ivory	3

For fitting tips and sizing chart see page 94 and 95.

*while stocks last until 2025

1) OVERVIEW OF CUP FORMS:

1 - Shallow cup form: concave silhouette, natural shape with additional volume in the lower cup area and a flatter upper breast; wider basic form

2 – Average cup form: even silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a larger underbust and small to medium cup size

3 – Full cup form: convex silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a smaller underbust and large cup size

2) OVERVIEW OF FORMS:

S = Symmetrical: Can be worn on the left and right side interchangeably



A = Asymmetrical: Specifically designed to be worn on the left or right side respectively. Extended form to compensate for missing tissue under the arm

E = Extra: Specifically designed to be worn on the left or right side respectively. Particularly suitable for women with a smaller cup size and fuller breast shape

*Light und Xtra Light forms have a more muted colour.

NATURA

Breast form with a temperature-regulating back.



NATURA double-layer breast form

-  Soft, matte PU film and premium-quality silicone guarantee a pleasant feel
-  **Comfort+**
Reduces heat and perspiration behind the breast form

PARTICULARLY SUITABLE FOR WOMEN WHO:

- have soft breast tissue (Natura 1S and 3E)
- want or need a very light breast form (Natura Xtra Light)

NATURA triple-layer breast form

-  Soft, matte PU film and premium-quality silicone guarantee a pleasant feel
-  **Comfort+**
Reduces heat and perspiration behind the breast form



Soft gel back

Sits comfortably against sensitive scar tissue

or

-  **Cosmetic**
Lightweight form that comes in plain, natural colours

PARTICULARLY SUITABLE FOR WOMEN WHO:

- have sensitive scar tissue

WEIGHT OPTIONS



COLOURS



Advanced

Overview

Image	Product name	Item no.	Cup form ¹	Weight	Shape ²	Sizes	Colour*	Layers
	Natura Xtra Light 1S ^N	401	1 – shallow	– 40%	S	1–14	Ivory	2
	Natura Xtra Light 2S ^N	400	2 – average	– 40%	S	1–17	Ivory Tawny	2
	Natura Light 1S ^N	402	1 – shallow	– 25%	S	1–14	Ivory	2
	Natura Light 1S	664	1 – shallow	– 25%	S	1–14	Ivory	3
	Natura Light 2S	390	2 – average	– 25%	S	1–17	Ivory Tawny	3
	Natura Light 2A	392	2 – average	– 25%	A	0–14 L/R	Ivory	3
	Natura Light 3S	391	3 – volume	– 25%	S	0–12	Ivory	3
	Natura Cosmetic 2S ^{N*}	323	2 – average	–20%	S	1–14	Ivory	3
	Natura Cosmetic 2S	320	2 – average	–20%	S	1–14	Ivory	3
	Natura Cosmetic 3S*	321	3 – volume	–20%	S	1–12	Ivory	3
	Natura Cosmetic 3E*	322	3 – volume	–20%	E	1–14 L/R	Ivory	3
	Natura 1S	396	1 – shallow	full even weight distribution	S	1–14	Ivory	2
	Natura 2S	398	2 – average	full even weight distribution	S	1–10	Ivory	3
	Natura 2A	393	2 – average	full even weight distribution	A	0–12 L/R	Ivory	3
	Natura 3S	395	3 – volume	full even weight distribution	S	0–12	Ivory	3
	Natura 3E	397	3 – volume	full even weight distribution	E	1–14 L/R	Ivory	2

For fitting tips and sizing chart see page 94 and 95.

*while stocks last until 2025

1) OVERVIEW OF CUP FORMS:

1 - Shallow cup form: concave silhouette, natural shape with additional volume in the lower cup area and a flatter upper breast; wider basic form

2 – Average cup form: even silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a larger underbust and small to medium cup size

3 – Full cup form: convex silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a smaller underbust and large cup size

2) OVERVIEW OF FORMS:

S = Symmetrical: Can be worn on the left and right side interchangeably


A = Asymmetrical: Specifically designed to be worn on the left or right side respectively. Extended form to compensate for missing tissue under the arm

E = Extra: Specifically designed to be worn on the left or right side respectively. Particularly suitable for women with a smaller cup size and fuller breast shape

*Light und Xtra Light forms have a more muted colour.

ESSENTIAL

Breast form that provides the basic care your patient needs for a symmetrical appearance



Single-layer breast form



Soft, matte PU film and high-quality silicone guarantee a pleasant feel

Standard and lightweight breast forms to create a symmetrical appearance and restore a feminine silhouette.

Basic

Overview

Image	Product name	Item no.	Cup form ¹	Weight	Shape ²	Sizes	Colour [*]	Layers
	Essential Light 1S^N	314	1 – shallow	– 25%	S	1–14	Ivory	1
	Essential Light 2S	442	2 – average	– 25%	S	1–17	Ivory Tawny	1
	Essential Light 3E	556	3 – volume	– 25%	E	1–14 L/R	Ivory Tawny (upon request)	1
	Essential 3S	363	3 – volume	full even weight distribution	S	0–12	Ivory	1

For fitting tips and sizing chart see page 94 and 95.

WEIGHT OPTIONS



Standard weight



-25%

COLOURS



Ivory



Tawny

1) OVERVIEW OF CUP FORMS:

1 - Shallow cup form: concave silhouette, natural shape with additional volume in the lower cup area and a flatter upper breast; wider basic form

2 – Average cup form: even silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a larger underbust and small to medium cup size

3 – Full cup form: convex silhouette, main volume concentrated in the upper and mid-cup area, especially suitable for women with a smaller underbust and large cup size

2) OVERVIEW OF FORMS:

S = Symmetrical: Can be worn on the left and right side interchangeably

A = Asymmetrical: Specifically designed to be worn on the left or right side respectively. Extended form to compensate for missing tissue under the arm

E = Extra: Specifically designed to be worn on the left or right side respectively. Particularly suitable for women with a smaller cup size and fuller breast shape

^{*}Light und Xtra Light forms have a more muted colour.