



SIGVARIS
GROUP

Medical Compression Therapy Solutions



2024

sigvaris

mobilis
by SIGVARIS GROUP

Our mission

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our mission for helping people feel their best is why we've been a global market leader for more than 155 years.

Our compression portfolio includes Medical, Consumer Health Care, products for Medically Complex Edema, and Soft Orthopedic Goods that are specifically designed to restore, assist and enhance people's lives. We are passionate about making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

WE'VE GOT YOU COVERED
FOR ALL YOUR COMPRESSION
THERAPY NEEDS!

Live your life. Be yourself.

We produce the best compression garments to enhance our customers' health and quality of life.

We are dedicated to helping legs look and feel their best.



RECIPROCATED TOE PRODUCTS

Crafted with a smooth, flat toe seam and contoured toe box for added comfort.



THE MOST COMFORTABLE FABRIC OPTIONS AVAILABLE

Including the finest cottons, wools, microfibers and double covered yarns.



WE ARE DEDICATED TO INNOVATION

We continue to invest in medical research to help people feel their best. Every day.



THE HIGHEST QUALITY STANDARDS

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality.



WE OFFER MORE VARIETY

More materials, colours, patterns, and styles than any other medical compression brand.



THE LARGEST SELECTION OF PREMIUM NATURAL FIBERS

Including Sea Island Cotton, Supima Cotton, Linen and Merino Wool.



DESIGNED FOR YOUR BEST FIT

Using studies on human data measurements, we fit 99.49% of the population.

WHY SIGVARIS GROUP?

| | | PRODUCT | STYLES | COMPRESSION LEVEL | PAGE |
|---------------------------------------|----------------------------------|--|---|--|---------------|
| MEDICAL COMPRESSION THERAPY SOLUTIONS | ESSENTIAL | COTTON 230 | <ul style="list-style-type: none"> ● calf, calf w/grip-top, thigh w/grip-top ○ calf, thigh w/grip-top | 20–30 mmHg 30–40 mmHg | 8 |
| | | OPAQUE 860 | <ul style="list-style-type: none"> ● calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose ○ calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh w/waist attachment | 20–30 mmHg 30–40 mmHg | 10, 12 |
| | | SHEER 780 <i>(for women only)</i> | <ul style="list-style-type: none"> ● calf, thigh w/grip-top, pantyhose, ○ calf, thigh w/grip-top, pantyhose | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 16 |
| | STYLE | MEDIUM SHEER 750 <i>(for women only)</i> | ● calf, thigh w/grip-top, pantyhose | 20–30 mmHg | 18 |
| | | SOFT OPAQUE 840 <i>(for women only)</i> | <ul style="list-style-type: none"> ● calf, calf w/grip-top, thigh w/grip-top, pantyhose, maternity pantyhose ○ calf, thigh w/grip-top, pantyhose | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 20 |
| | | SEA ISLAND COTTON 220 | ● calf | 20–30 mmHg | 22 |
| | | MERINO WOOL 240 | ● calf | 20–30 mmHg | 24 |
| | | LINEN 250 | ● calf | 20–30 mmHg | 26 |
| | | MICROFIBER 820 <i>(for men only)</i> | ● calf, calf w/grip-top, thigh w/grip-top (extended foot sizes available) | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 28 |
| | | MICROFIBER PATTERNS 830 | ● calf | 20–30 mmHg | 30 |
| | | MOTION | CUSHIONED COTTON 360 | ● calf | 20–30 mmHg |
| | HIGH TECH 412C | | ● calf | 20–30 mmHg | 36 |
| | THERMOREGULATING WOOL 422 | | ● calf | 20–30 mmHg | 38 |
| | SPECIALTY | NATURAL RUBBER 500 | ○ calf, thigh, thigh w/grip-top, thigh w/waist attachment, pantyhose | 30–40 mmHg 40–50 mmHg 50–60 mmHg | 42 |
| | | ADVANCE ARMSLEEVE 910 | ○ armsleeve w/grip-top w/gauntlet and without gauntlet | 15–20 mmHg 20–30 mmHg | 44 |
| | INELASTIC WRAPS | | | 46 | |

Closed Toe = ● Open Toe = ○

| | PRODUCT | STYLES | COMPRESSION LEVEL | PAGE |
|------------------------------------|---|---|-------------------|---------------|
| CONSUMER HEALTH | CASUAL COTTON 146/186 | ● calf | 15–20 mmHg | 66 |
| | TRAVENO TRAVEL SOCKS 146/T - 186/T | ● calf | 15–20 mmHg | 67 |
| | SEA ISLAND COTTON 151/191 | ● calf | 15–20 mmHg | 68 |
| | ALL-SEASON MERINO WOOL 152/192 | ● calf | 15–20 mmHg | 70 |
| | LINEN 153/193 | ● calf | 15–20 mmHg | 71 |
| | MICROFIBER SHADES 143/183 | ● calf | 15–20 mmHg | 72 |
| | SHEER FASHION 120 | ● calf, thigh w/grip-top, pantyhose ○ calf | 15–20 mmHg | 74 |
| | SOFT SILHOUETTE LEGGINGS 170 | ○ open foot | 15–20 mmHg | 76 |
| | ATHLETIC RECOVERY SOCKS 401 | ● calf | 15–20 mmHg | 78 |
| | PERFORMANCE SLEEVES 412V | ○ open foot | 20–30 mmHg | 79 |
| | MERINO OUTDOOR SOCKS 421 | ● calf | 15–20 mmHg | 80 |
| EVERSOFT DIABETIC SOCKS 160 | ● calf | 8–15 mmHg | 81 | |
| | ACCESSORIES | | | 82 |
| RESOURCES | SYMPTOMS & INDICATIONS | | | 90 |
| | GLOSSARIES OF ICONS & TERMS | | | 91-92 |
| | HOW TO MEASURE & ORDER | | | 93-95 |
| MOBILIS | mobilis ORTHOPEDIC BRACING | | | 96-117 |
| | ELBOW | | | 100 |
| | WRIST & THUMB | | | 103 |
| | BACK | | | 107 |
| | KNEE | | | 110 |
| | ANKLE | | | 114 |
| | HOW TO ORDER MOBILIS | | | 117 |

sigvaris ESSENTIAL

Medical Compression Therapy Solutions



Support Yourself

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

ESSENTIAL PRODUCTS

| Series | Product | Page |
|------------|----------------|------|
| 230 | Cotton | 8 |
| 860 | Opaque (Women) | 10 |
| 860 | Opaque (Men) | 12 |

BEST SUITED FOR

- Containment of various types of edema, including venous edema, phlebolymphe~~ma~~, and early stages of lymphedema
- Ideal for dress or casual shoes
- Suited for post-vein treatment

Featured colour is Light Beige

COTTON 230

for women & men

Compression levels

20–30 mmHg

30–40 mmHg

Features and Benefits



Sensitive Skin



Cotton



Breathable

- 100% Supima Cotton against the skin
- Ideal for sensitive skin and allergies
- Breathable and durable
- Excellent containment of various types of edema
- Classic ribbed design for everyday wear

Fiber Content

- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Calf w/ Grip-top
Medium Long in Black

232CMLM99/S



Silicone dotted grip-top available for calf and thigh

SIZING CHART

| | Small | | Medium | | Large | | X Large | |
|---|-------------------------|-------------------|-----------------------------|-------------------|---------------------------|-------------------|-------------------------|-----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | | NA | |
| Leg Length: Women & Open Toe | Short | Long | Short | Long | Short | Long | | |
| | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | | |
| Leg Length: Men | Short | Long | Short | Long | Short | Long | | |
| | ≤76 cm (30") | >76 cm (30") | ≤76 cm (30") | >76 cm (30") | ≤76 cm (30") | >76 cm (30") | | |
| Size | SS | SL | MS | ML | LS | LL | XS | XL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

Open Toe = ○



| | | | | | |
|---|---|---|---|-----|---|
| Calf (C)* | ● | ● | ● | ●● | ● |
| Calf w/Grip-Top (C/S) Silicone dotted grip-top band | | | | ●** | ● |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | | | | ●● | ● |

* 30–40 mmHg extra large calf only available in black with grip-top

** 20–30 mmHg only

Product Code Example: 232NMLW66

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 23 | 2 | N | ML | W | 66 |

Product Code Example: 232CLSM66/S

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour | Grip-top |
|--------|------------|-------|------|-----------------|--------|----------|
| 23 | 2 | C | LS | M | 66 | /S |

OPAQUE 860

for women

Compression Levels

20–30 mmHg

30–40 mmHg

Features and Benefits



Durable



Largest Size Range



Opaque

- Largest size range, from petite to plus
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema
- Plus size pantyhose are an alternative to maternity pantyhose

Fiber Content

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free

Featured colour is Honey



20–30 mmHg
Closed Toe Calf Medium Long
in Black

862CMLW99



Silicone dotted grip-top available
for calf and thigh

SIZING CHART

| | Small | | | Medium | | | Large | | X Large | |
|---|---------------------------|----------------|----------------|--------------------------|------------------------|----------------|---------------------------|----------------|----------------------|----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | | 21.5–25.5 cm (8.5"–10") | | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | | 31.5–45 cm (12.5"–17.5") | | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Circumference: Calf Plus w/Grip-Top | NA | | | NA | 39.5–53 cm (15.5"–21") | | 45.5–61 cm (18"–24") | | 51–71 cm (20"–28") | |
| Calf Length | Petite | Short | Long | Petite | Short | Long | Short | Long | Short | Long |
| | ≤35.5 cm (14") | ≤41 cm (16") | >41 cm (16") | ≤35.5 cm (14") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | | 48–73 cm (19"–28.5") | | | 52–84 cm (20.5"–33") | | 64–95 cm (25"–37.5") | |
| Leg Length | Petite | Short | Long | Petite | Short | Long | Short | Long | Short | Long |
| | ≤65 cm (25.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤65 cm (25.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") |
| Hip Circumference | 83–112.5 cm (32.5"–44.5") | | | 87–127 cm (34.5"–50") | | | 98–150.5 cm (38.5"–59.5") | | NA | |
| Plus or Maternity Size Hip Circumference | NA | ≤142 cm (56") | | NA | ≤157.5 cm (62") | | ≤178 cm (70") | | NA | |
| Size | SP | SS | SL | MP | MS | ML | LS | LL | XS | XL |
| Calf Plus w/Grip-Top | | | | | MS/P | ML/P | LS/P | LL/P | XS/P | XL/P |

Women's Shoe Size : Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

| HONEY 33 | GOLDEN 36 | LIGHT BEIGE 66 | BLACK 99 |
|-------------|--------------|----------------------|-------------|
|-------------|--------------|----------------------|-------------|

| | | | | |
|---|---|---|----|-----|
| Calf (C) | ● | ● | ○ | ●● |
| Calf Petite Sizes (C)* Not available with dotted grip-top band | ● | ● | | ● |
| Calf w/Grip-Top (C/S) Silicone dotted grip-top band | | | ●● | ● |
| Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band | | | ●● | ●●* |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | ● | ● | ○ | ●● |
| Thigh w/Grip-Top (N) XS & XL Silicone dotted grip-top band | | | ●● | ○ |
| Thigh Petite Sizes (N)* | ● | ● | | ● |
| Pantyhose (P)** | ● | ● | | ● |
| Pantyhose Petite (P)* | ● | ● | | ● |
| Pantyhose Plus and Maternity (P/P)** | | | ● | ● |
| Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R) | | | ○ | |

* 20-30 mmHg only

** Not available in XS or XL sizing

Product Code Example: 863NMSW33

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour |
|--------|------------|-------|------|------------------|--------|
| 86 | 3 | N | MS | W | 33 |

Product Code Example: 862CLLO66/P

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Plus size with Grip-top |
|--------|------------|-------|------|------------------|--------|-------------------------|
| 86 | 2 | C | LL | O | 66 | /P |



Featured colour is Ribbed Black

OPAQUE 860

for men

Compression Levels

20–30 mmHg / 30–40 mmHg

Features and Benefits



Durable



Largest Size Range



Opaque

- Durable construction with a ribbed design on closed toe calf products
- Smooth opaque design on open-toe styles
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema

Fiber Content

- 64% Nylon, 36% Spandex (content can vary slightly by style)
- Latex free



20–30 mmHg
Open Toe Thigh Medium Long
in Light Beige

862NMLO66



Silicone dotted grip-top available for calf and thigh

SIZING CHART

| | Small | | Medium | | Large | | X Large | |
|---|-------------------------|-------------------|-----------------------------|-------------------|---------------------------|-------------------|-------------------------|-------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Circumference: Calf Plus w/Grip-Top | NA | | 39.5–53 cm (15.5"–21") | | 45.5–61 cm (18"–24") | | 51–71 cm (20"–28") | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | | 64–95 cm (25"–37.5") | |
| Leg Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") |
| Size | SS | SL | MS | ML | LS | LL | XS | XL |
| Calf Plus w/Grip-Top | | | MS/P | ML/P | LS/P | LL/P | XS/P | XL/P |

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○



| | | | | |
|--|----|----|---|----|
| Calf (C) | ●* | ●* | ○ | ○ |
| Calf w/Grip-Top (C/S) Silicone dotted grip-top band | | | ○ | |
| Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band | | | ○ | ○* |
| Open Toe Thigh w/Grip-Top (N) Silicone dotted grip-top band | | | ○ | ○ |
| Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R) | | | ○ | |

* 20-30 mmHg only

** Not available in XS or XL with waist attachment

Product Code Example: 862CMLM99

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 86 | 2 | C | ML | M | 99 |

Product Code Example: 862NMLO66

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 86 | 2 | N | ML | O | 66 |

Medical Compression Therapy Solutions



Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

STYLE PRODUCTS

| Series | Product | Page |
|------------|---------------------|------|
| 780 | Sheer | 16 |
| 750 | Medium Sheer | 18 |
| 840 | Soft Opaque | 20 |
| 220 | Sea Island Cotton | 22 |
| 240 | Merino Wool | 24 |
| 250 | Linen | 26 |
| 820 | Microfiber | 28 |
| 830 | Microfiber Patterns | 30 |

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress attire



Featured colour is Toasted Almond

SHEER 780

for women

Compression Levels



Features and Benefits



Sheer



Elegant



Breathable

- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Excellent breathability and moisture management

Fiber Content

- 15-20 and 20-30 mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30-40 mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free



15-20 mmHg
Closed Toe Thigh Medium Long
in Dark Navy

781NMLW08



Sensinnov® grip-top on thigh

SIZING CHART

| | Small | | Medium | | Large | |
|----------------------------|------------------------------|-------------------|-----------------------------|-------------------|------------------------------|-------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| Calf Length | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | |
| Leg Length | Short | Long | Short | Long | Short | Long |
| | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") |
| Hip Circumference | 83–112.5 cm (32.5"–44.5") | | 87–127 cm (34.5"–50") | | 98–150.5 cm (38.5"–59.5") | |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

| | DARK NAVY 08 | HONEY 33 | WARM SAND 34* | GOLDEN 36 | TOASTED ALMOND 37* | CAFÉ 73* | MOCHA 85* | NIGHT-SHADE 94 | BLACK 99 |
|--|--------------|----------|---------------|-----------|--------------------|----------|-----------|----------------|----------|
| Calf (C) | ● | ●● | ●● | ●● | ●● | ●● | ●● | ● | ●● |
| Thigh w/Grip-Top (N) Sensinnov® grip-top band | ● | ●● | ●● | ●● | ●● | ●● | ●● | ● | ●● |
| Pantyhose (P) | ● | ●● | ●● | ●● | ●● | ●● | ●● | ● | ●● |

* 20–30 mmHg only

Product Code Example: 781NMLW36

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 78 | 1 | N | ML | W | 36 |

Product Code Example: 783PSSO33

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 78 | 3 | P | SS | O | 33 |



Featured colour is Black

MEDIUM SHEER 750

for women

Compression Level

20-30 mmHg

Features and Benefits



Medium Sheer



Durable



Breathable

- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Comfortable 2" top band on calf style

Fiber Content

- 68% Nylon, 32% Spandex (content can vary slightly by style)
- Latex free



20-30 mmHg
Closed Toe Thigh Medium Long in Honey

752NMLW33



Silicone dotted grip-top on thigh

SIZING CHART

| | Small | | Medium | | Large | |
|----------------------------|------------------------------|-------------------|-----------------------------|-------------------|------------------------------|-------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| Calf Length | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | |
| Leg Length | Short | Long | Short | Long | Short | Long |
| | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") |
| Hip Circumference | 83–112.5 cm (32.5"–44.5") | | 87–127 cm (34.5"–50") | | 98–150.5 cm (38.5"–59.5") | |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

| | | | |
|--|---|---|--|
| Closed Toe = ● |  |  |  |
| | HONEY 33 | NIGHT-SHADE 94 | BLACK 99 |
| Calf (C) | ● | ● | ● |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | ● | ● | ● |
| Pantyhose (P) | ● | ● | ● |

Product Code Example: 752NMLW33

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 75 | 2 | N | ML | W | 33 |



Featured colour is Mulberry

SOFT OPAQUE 840

for women

Compression Levels



Features and Benefits



Opaque



Microfiber



Soft

- Fashionable opaque hosiery for professional and everyday wear
- Long-lasting microfiber softness provides all day comfort and allows product to easily glide on and off
- Durable and breathable
- Maternity pantyhose option available

Fiber Content

- 79% Nylon, 21% Spandex (content can vary slightly by style)
- Latex free



20-30 mmHg
Closed Toe Pantyhose Medium Long in Black

842PMLW99



Silicone dotted grip-top for calf and Sensinno[®] grip-top for thigh



SIZING CHART

| | Small | | Medium | | Large | |
|------------------------------------|------------------------------|-------------------|-----------------------------|-------------------|------------------------------|-------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| Calf Length | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | |
| Leg Length | Short | Long | Short | Long | Short | Long |
| | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") |
| Hip Circumference | 83–112.5 cm (32.5"–44.5") | | 87–127 cm (34.5"–50") | | 98–150.5 cm (38.5"–59.5") | |
| Maternity hip circumference | ≤142 cm (56") | | ≤157.5 cm (62") | | ≤178 cm (70") | |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

| | MIDNIGHT BLUE 09 | CHAI 35 | PECAN 40 | MULBERRY 60* | ESPRESSO 89* | GRAPHITE 91 | BLACK 99 |
|--|------------------------|------------|-------------|-----------------|-----------------|----------------|-------------|
|--|------------------------|------------|-------------|-----------------|-----------------|----------------|-------------|

| | | | | | | | |
|---|---|----|---|---|---|---|----|
| Calf (C) | ● | ●○ | ● | ● | ● | ● | ●○ |
| Calf w/Grip-Top (C/S)* Silicone dotted grip-top | | | | | | | ● |
| Thigh w/Grip-Top (N) Sensinnov® grip-top band | ● | ●○ | ● | ● | ● | ● | ●○ |
| Pantyhose (P) | ● | ●○ | ● | ● | ● | ● | ●○ |
| Maternity Pantyhose (M)* | | | | | | | ● |

* 20–30 mmHg only

Product Code Example: 842CMSW60

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 84 | 2 | C | MS | W | 60 |



SEA ISLAND COTTON 220

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 68

Features and Benefits



Soft



Cotton



Breathable

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free

Featured colour is Black



20–30 mmHg
Closed Toe Calf Small Long
in Brown

222CSLW11



20–30 mmHg
Closed Toe Calf Medium Long
in Navy

222CMLM10

SIZING CHART

| | Small | | Medium | | Large | | X Large | |
|----------------------------|-------------------------|-----------------|-----------------------------|-----------------|---------------------------|-----------------|-------------------------|-----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL | XS | XL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



| | | | |
|-----------------|---|---|---|
| Calf (C) | ● | ● | ● |
|-----------------|---|---|---|

Product Code Example: 222CLLM10

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 22 | 2 | C | LL | M | 10 |

MERINO WOOL 240

for women & men

Compression Level

20-30 mmHg

Also available in 15-20 mmHg, see page 70

Features and Benefits



Wool



Thermoregulating



Stylish



- High quality extra-fine Merino wool makes this sock soft and breathable
- Moisture wicking capabilities keep feet dry
- Naturally odour resistant
- Thermoregulating for year-round wear
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



Featured colour is Black



20-30 mmHg
Closed Toe Calf Medium Long
in Brown

242CMLM11



20-30 mmHg
Closed Toe Calf Medium Long
in Navy

242CMLW10

SIZING CHART

| | Small | | Medium | | Large | | X Large | |
|----------------------------|-------------------------|-----------------|-----------------------------|-----------------|---------------------------|-----------------|-------------------------|-----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL | XS | XL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



| | | | |
|-----------------|---|---|---|
| Calf (C) | ● | ● | ● |
|-----------------|---|---|---|

Product Code Example: 242CMSW99

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 24 | 2 | C | MS | W | 99 |

Featured
colour is
Denim

LINEN 250

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 71

Features and Benefits



Hypoallergenic
Natural fibers



Anti-bacterial



Moisture
Management

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Long
in Light Grey

252CMLM205



20–30 mmHg
Closed Toe Calf Medium Long
in Lavender

252CMLM207

SIZING CHART

| | Small | | Medium | | Large | |
|----------------------------|-------------------------|-----------------|-----------------------------|-----------------|---------------------------|-----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| Calf Length | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



| | | | | |
|-----------------|---|---|---|---|
| Calf (C) | ● | ● | ● | ● |
|-----------------|---|---|---|---|

Product Code Example: 252CMLM207

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 25 | 2 | C | ML | M | 207 |

MICROFIBER 820

for men

Compression Levels

15-20 mmHg

20-30 mmHg

30-40 mmHg

Features and Benefits



Microfiber



Soft



Stylish

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking
- Glides on and off with ease
- Comfortable and durable sock for professional and everyday wear
- Available in extended foot sizes
- Smooth, flat toe seam

Fiber Content


- 87% Nylon , 13% Spandex
(content can vary slightly by style)
- Latex free

Featured colour is Tan Khaki



20-30 mmHg
Closed Toe Thigh Medium Long in Black

822NMLM99



Silicone dotted grip-top option available for calf. Sensinno[®] grip-top for thigh.



SIZING CHART

| | Small | | Medium | | Large | | X Large | |
|----------------------------|-------------------------|-----------------|-----------------------------|------------------|---------------------------|-----------------|-------------------------|------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | | NA | |
| Leg Length | Short | Long | Short | Long | Short | Long | | |
| | ≤76 cm (30") | >76 cm (30") | ≤76 cm (30") | >76 cm (30") | ≤76 cm (30") | >76 cm (30") | | |
| Size Shoe Size | SS 7-11 | SL 9.5-14 | MS 7-11 | ML 9.5-14 | LS 7-11 | LL 9.5-14 | XS 7-11 | XL 9.5-14 |
| | | | | MX* 12.5-17.5 | | | LX* 12.5-17.5 | XX* 12.5-17.5 |

Men's shoe sizes: Short Leg 7–11, Long Leg 9.5–14 (For extended foot sizes: long leg 12.5–17.5)

*MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30 mmHg only. Extended foot sizes not available with grip top

STYLES AND COLOURS

Closed Toe = ●



| | TAN KHAKI 32 | STEEL GREY 93 | BLACK 99 |
|---|-----------------|------------------|-------------|
| Calf (C) | ● | ●*** | ● |
| Calf w/Grip-Top (C/S)** Silicone dotted grip-top band | ● | | ● |
| Thigh w/Grip-Top (N)** Sensinnov® grip-top band | ● | | ● |

** Not available in 15-20 mmHg

*** Steel Grey only available in 15-20 mmHg and 20-30 mmHg

Product Code Example: 823CMLM32

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 82 | 3 | C | ML | M | 32 |

Product Code Example: 822CMLM32/S

| Series | mmHg Level | Style | Size | Gender | Colour | Grip-top |
|--------|------------|-------|------|--------|--------|----------|
| 82 | 2 | C | ML | M | 32 | /S |

MICROFIBER PATTERNS 830

for women & men

Compression Level

20-30 mmHg

Also available in 15-20 mmHg, see page 72

Features and Benefits



Microfiber



Stylish



Patterns

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Fashionable, durable sock for everyday wear
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex (content can vary slightly by style)
- Latex free



Featured colour is Wisteria Windows



20-30 mmHg
Closed Toe Calf Small Long
in Royal Blue Argyle

832CSLM48



20-30 mmHg
Closed Toe Calf Small Short
in Mariner Stripe

832CSSW02

SIZING CHART

| | Small | | Medium | | Large | |
|----------------------------|-------------------------|-----------------|-----------------------------|-----------------|---------------------------|-----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| Calf Length | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12
 Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Women

Closed Toe = ●



| | | | | | |
|-----------------|---|---|---|---|---|
| Calf (C) | ● | ● | ● | ● | ● |
|-----------------|---|---|---|---|---|



| | | | | | |
|-----------------|---|---|---|---|---|
| Calf (C) | ● | ● | ● | ● | ● |
|-----------------|---|---|---|---|---|

Men

Closed Toe = ●



| | | | | | |
|-----------------|---|---|---|---|---|
| Calf (C) | ● | ● | ● | ● | ● |
|-----------------|---|---|---|---|---|



| | | |
|-----------------|---|---|
| Calf (C) | ● | ● |
|-----------------|---|---|

Product Code Example: 832CMLW02

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 83 | 2 | C | ML | W | 02 |

sigvaris MOTION

Medical Compression Therapy Solutions



Push Yourself

Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool (422)
- Moisture-wicking yarns

MOTION PRODUCTS

| Series | Product | Page |
|-------------|-----------------------|------|
| 360 | Cushioned Cotton | 34 |
| 412C | High Tech | 36 |
| 422 | Thermoregulating Wool | 38 |

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Casual and occupational footwear
- Walking and everyday activities

CUSHIONED COTTON 360

for women & men

Compression Level

20-30 mmHg

Features and Benefits



Cushioned Sole



Cotton



Sensitive Skin

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Smooth, flat toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time

Fiber Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free

Featured colour is White



20-30 mmHg
Closed Toe Calf Medium Long
in Black

362CMLW99



Flat, low friction toe seam

SIZING CHART

| | Small | | Medium | | Large | | X Large | |
|----------------------------|-------------------------|-----------------|-----------------------------|-----------------|---------------------------|-----------------|-------------------------|-----------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL | XS | XL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



| | | |
|-----------------|---|---|
| Calf (C) | ● | ● |
|-----------------|---|---|

Product Code Example: 362CMLM00

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|------------|-------|------|--------|--------|
| 36 | 2 | C | ML | M | 00 |



Featured colour is Lime

HIGH TECH 412C

for women & men

Compression Level

20-30 mmHg

Features and Benefits



Moisture Wicking



Anatomically Contoured



Colours

- High-Tech fibers for moisture wicking and improved thermoregulation
- Cushioned footbed, toe, and Achilles tendon zones for optimal comfort
- Anatomically contoured to support the natural shape of the foot and calf
- Smooth, flat toe seam
- Perfect for everyday activities and active individuals

Fiber Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free



20-30 mmHg
Closed Toe Calf Small Medium in Steel Blue

412CSM62



20-30 mmHg
Closed Toe Calf Medium Long in Pink


412CML56

SIZING CHART

| | Small | | | Medium | | | Large | | | X Large | | |
|----------------------------|-------------------------|--------|---------|-----------------------------|--------|---------|---------------------------|---------|---------|-------------------------|---------|---------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | | 21.5–25.5 cm (8.5"–10") | | | 25.5–29 cm (10"–11.5") | | | 29–33 cm (11.5"–13") | | |
| Calf Circumference | 28–39 cm (11"–15.5") | | | 31.5–45 cm (12.5"–17.5") | | | 35.5–51 cm (14"–20") | | | 40.5–58 cm (16"–23") | | |
| Women's Shoe Size | 4.5–7.5 | 8–10.5 | 11–13.5 | 4.5–7.5 | 8–10.5 | 11–13.5 | 8–10.5 | 11–13.5 | >14 | 8–10.5 | 11–13.5 | >14 |
| Men's Shoe Size | ≤6.5 | 7–9.5 | 10–12.5 | ≤6.5 | 7–9.5 | 10–12.5 | 7–9.5 | 10–12.5 | 13–15.5 | 7–9.5 | 10–12.5 | 13–15.5 |
| Size | SS | SM | SL | MS | MM | ML | LS | LM | LL | XS | XM | XL |

STYLES AND COLOURS

Closed Toe = ●



| | | | | | | | |
|-----------------|---|---|---|---|---|---|---|
| Calf (C) | ● | ● | ● | ● | ● | ● | ● |
|-----------------|---|---|---|---|---|---|---|

Product Code Example: 412CML54

| Series | mmHg Level | Style | Size | Colour |
|--------|------------|-------|------|--------|
| 41 | 2 | C | ML | 54 |



Featured colour is Charcoal

THERMO-REGULATING WOOL 422

for women & men

Compression Level

20-30 mmHg

Also available in 15-20 mmHg, see page 80

Features and Benefits



Wool



Cushioned Sole



Thermoregulating

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned, padded soles help prevent blisters
- Soft and comfortable
- Smooth, flat toe seam
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free



20-30 mmHg
Closed Toe Calf Large Long
in Charcoal

422CLL12



20-30 mmHg
Closed Toe Calf Medium Short
in Olive

422CMS31

SIZING CHART

| | Small | | | Medium | | | Large | | |
|----------------------------|-------------------------|--------|---------|-----------------------------|--------|---------|---------------------------|---------|---------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | | 21.5–25.5 cm (8.5"–10") | | | 25.5–29 cm (10"–11.5") | | |
| Calf Circumference | 28–39 cm (11"–15.5") | | | 31.5–45 cm (12.5"–17.5") | | | 35.5–51 cm (14"–20") | | |
| Women's Shoe Size | 4.5–7.5 | 8–10.5 | 11–13.5 | 4.5–7.5 | 8–10.5 | 11–13.5 | 8–10.5 | 11–13.5 | >14 |
| Men's Shoe Size | ≤6.5 | 7–9.5 | 10–12.5 | ≤6.5 | 7–9.5 | 10–12.5 | 7–9.5 | 10–12.5 | 13–15.5 |
| Size | SS | SM | SL | MS | MM | ML | LS | LM | LL |

STYLES AND COLOURS

Closed Toe = ●



| | | |
|-----------------|---|---|
| Calf (C) | ● | ● |
|-----------------|---|---|

Product Code Example: 422CML12

| Series | mmHg Level | Style | Size | Colour |
|--------|------------|-------|------|--------|
| 42 | 2 | C | ML | 12 |

sigvaris SPECIALTY

Medical Compression Therapy Solutions



Take care of Yourself

Why choose **Specialty**?

- Specific medical needs
- Need low fatigue fibers
- Natural rubber (500)

SPECIALTY PRODUCTS

| Series | Product | Page |
|------------|-------------------|------|
| 500 | Natural Rubber | 42 |
| 910 | Advance Armsleeve | 44 |

BEST SUITED FOR

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Early stages of lymphedema
- Active venous ulcers and prevention of ulcer recurrence
- People with early stages of arm lymphedema
- Post-mastectomy care



NATURAL RUBBER 500

for women & men

Compression Levels

30–40 mmHg

40–50 mmHg

50–60 mmHg

Features and Benefits



Natural Rubber



Durable



Largest Size Range

- Durable, long-term containment solution with minimal fiber fatigue throughout the day
- Extended sizing for a custom-like fit
- **Ideal for patients with more serious indications, including:**
 - Severe varicosities
 - Chronic edema post DVT
 - Prevention and treatment of Active Venous Ulcers
 - Severe venous stasis
 - Lymphedema

Fiber Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction



30–40 mmHg
Open Toe Thigh Medium Average
Long in Beige

503NM2077



Silicone dotted grip-top available
on thigh

SIZING CHART

| | Small | | | | Medium | | | | Large | | | | X Large | | | |
|---|-----------------------|-----------------|----------------------|-----------------|-----------------------|-----------------|----------------------|-----------------|-----------------------|-----------------|----------------------|-----------------|---------------------------|-----------------|-----------------|-----------------|
| Ankle Circumference | 19–23 cm (7.5"–9") | | | | 23–30 cm (9"–12") | | | | 30–40 cm (12"–15.5") | | | | 40–46 cm (15.5"–18") | | | |
| Calf Circumference | Average | | Full | | Average | | Full | | Average | | Full | | Average | | Full | |
| | 28–38 cm (11"–15") | | ≤43 cm (17") | | 32–43 cm (12.5"–17") | | ≤48 cm (19") | | 36–48 cm (14.5"–19") | | ≤57 cm (22.5") | | 45–57 cm (18"–22.5") | | ≤68 cm (27") | |
| Calf Circumference for 505C (Full only) | 505C: full calf only | | 32–43 cm (12.5"–17") | | 505C: full calf only | | 36–48 cm (14.5"–19") | | 505C: full calf only | | 45–57 cm (18"–22.5") | | 505C: not available in XL | | | |
| Calf Length | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long |
| | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | >40 cm (>15.5") |
| Thigh Circumference | Average | | Full | | Average | | Full | | Average | | Full | | N/A | | | |
| | 41–54 cm (16"–21") | | ≤61 cm (24") | | 48–61 cm (19"–24") | | ≤71 cm (28") | | 58–71 cm (23"–28") | | ≤81 cm (32") | | | | | |
| Thigh Circumference for 505T (Full only) | 505T: full thigh only | | ≤61 cm (24") | | 505T: full thigh only | | ≤71 cm (28") | | 505T: full thigh only | | ≤81 cm (32") | | | | | |
| Leg Length | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | | | | |
| | ≤71 cm (28") | >71 cm (28") | ≤71 cm (28") | >71 cm (28") | ≤71 cm (28") | >71 cm (28") | ≤71 cm (28") | >71 cm (28") | ≤71 cm (28") | >71 cm (28") | ≤71 cm (28") | >71 cm (28") | | | | |
| Hip Circumference | <122 cm (48") | | | | <137 cm (54") | | | | <155 cm (61") | | | | | | | |
| Size | S1 | S2 | S3 | S4 | M1 | M2 | M3 | M4 | L1 | L2 | L3 | L4 | X1 | X2 | X3 | X4 |

STYLES AND COLOURS

Open Toe = ○



Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 90 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.

| | | |
|---|---|---|
| Calf (C) | ○ | Calf only available in S, M, and L full circumference for 50–60 mmHg |
| Thigh w/out Grip-Top (T) | ○ | |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | ○ | Not available in 50-60 mmHg or XL |
| Pantyhose (P) | ○ | Pantyhose only available in 30–40 mmHg and only in full thigh circumference |
| Thigh w/Waist Attachment (W) Indicate Leg LEFT (/L) or RIGHT (/R) | ○ | Thigh w/ waist attachment is available in average and full thigh circumferences in 30–40 mmHg and in full thigh circumferences only in 40–50 mmHg |

Product Code Example: 504CS2077

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour |
|--------|------------|-------|------|-----------------|--------|
| 50 | 4 | C | S2 | O | 77 |

Product Code Example: 503WL1077/L

| Series | mmHg Level | Style | Size | Gender/Open Toe | Colour | Leg |
|--------|------------|-------|------|-----------------|--------|-----|
| 50 | 3 | W | L1 | O | 77 | /L |



ADVANCE ARMSLEEVE 910

for women & men

Compression Levels

15-20 mmHg

20-30 mmHg

Features and Benefits



Moisture
Wicking



SensinnoV
Grip-top Band




Breathable

- Ideal for patients with arm lymphedema or post-mastectomy care
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable and soft
- SensinnoV® grip-top band ensures armsleeve stays in place all day long
- Available with or without gauntlet

Fiber Content

- 50% Nylon, 28% Rayon, 22% Spandex
- Latex free



With
Gauntlet



SensinnoV® grip-top band

20-30 mmHg
Armsleeve With Gauntlet
Small Long in Beige

912ASLG77/S

SIZING CHART

| | Small | | Medium | | | | Large | | | | X Large | | | |
|--------------------------|-------------------------|-----------------|---------------------------|-----------------|-------------------------|-----------------|---------------------------|-----------------|-------------------------|-----------------|---------------------------|-----------------|-------------------------|-----------------|
| Palm (A) | 17-21 cm (6.5"-8") | | 18-22 cm (7"-8.5") | | | | 19-23 cm (7.5"-9") | | | | 20-24 cm (8"-9.5") | | | |
| Wrist (C) | 14-16 cm (5.5"-6.5") | | 16-18 cm (6.5"-7") | | | | 18-20 cm (7"-8") | | | | 20-23 cm (8"-9") | | | |
| Upper Mid-Arm (F) | Regular | | Regular | | Plus | | Regular | | Plus | | Regular | | Plus | |
| | 26-34 cm (10"-13.5") | | 29-37 cm (11.5"-14.5") | | 37-43 cm (14.5"-17") | | 32-40 cm (12.5"-15.5") | | 40-46 cm (15.5"-18") | | 34-42 cm (13.5"-16.5") | | 42-48 cm (16.5"-19") | |
| Arm Length (C-G)* | Regular | Long | Regular | Long | Regular | Long | Regular | Long | Regular | Long | Regular | Long | Regular | Long |
| | ≤43 cm (17") | >43 cm (17") | ≤43 cm (17") | >43 cm (17") | ≤43 cm (17") | >43 cm (17") | ≤43 cm (17") | >43 cm (17") | ≤43 cm (17") | >43 cm (17") | ≤43 cm (17") | >43 cm (17") | ≤43 cm (17") | >43 cm (17") |
| Size | SR | SL | MR | ML | MR/P | ML/P | LR | LL | LR/P | LL/P | XR | XL | XR/P | XL/P |

* Should have **slight bend** in arm when measuring length


STYLES AND COLOURS

Open Sleeve = ○



| | |
|--|---|
| Armsleeve w/Gauntlet (G) Sensinnov® grip-top band | ○ |
| Armsleeve w/out Gauntlet (O) Sensinnov® grip-top band | ○ |

Note: all armsleeves come with Grip-top (/S)
15-20 mmHg not available with gauntlet or plus sizing



ARMSLEEVE WITHOUT GAUNTLET

Product Code Example: 912ASLO77/S

| Series | mmHg Level | Style | Size | W/out Gauntlet | Colour | Grip-Top |
|--------|------------|-------|------|----------------|--------|----------|
| 91 | 2 | A | SL | O | 77 | /S |

ARMSLEEVE WITH GAUNTLET

Product Code Example: 912ALRG77/PS

| Series | mmHg Level | Style | Size | W/Gauntlet | Colour | Plus Size | Grip-Top |
|--------|------------|-------|------|------------|--------|-----------|----------|
| 91 | 2 | A | LR | G | 77 | /P | S |

sigvaris INELASTIC WRAPS

Ready-to-wear wraps for the active treatment and management of conditions impacting the lymphatic and venous systems



Why Choose Inelastic compression wraps?

- Adjustable pressure
- Large variety of solutions

MEDICALLY COMPLEX EDEMA PRODUCTS

| Product | Page |
|----------------------------------|------|
| Coolflex™ | 48 |
| Compreflex™ | 52 |
| Compreflex™ Standard Calf & Foot | 54 |
| Compreflex™ Liner | 56 |
| Compreboot® | 57 |
| Compreflex™ Strap Extender | 58 |
| Accutab® | 58 |
| Compreflex™ Thigh | 59 |
| Compreflex™ Knee | 60 |
| Compreknee™ | 61 |
| Compreshort® | 62 |

The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems.

The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous disorders in the lower extremities. Sigvaris is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.

COOLFLEX™

Wrap Details

LIGHTWEIGHT, BREATHABLE COMPRESSION WEAR

Coolflex™ is a lightweight compression device, optimal for individuals with fragile or broken skin. Made of a non-stretch mesh material, the Coolflex™ transfers heat and wicks away moisture, leaving the skin cool and dry.

The features of the closure system, along with the built-in guidelines for application, promote independence and encourage improved patient outcomes.

BREATHABLE GARMENT

The Spacer Fabric® TR is designed to be highly breathable for cool-wearing comfort

ONE-HANDED DONNING

EZ-Close Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

INNOVATIVE DESIGN

Durable, metal hook closures for maximum strength and durability



20–50 mmHg



Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Well-managed Lymphedema
- Non-specific or post-surgical edema
- Easy to put on for patients with limited dexterity


Features

- Lightweight and breathable for cool-wearing comfort
- Durable metal clasps, instead of Velcro closure
- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort



ACCURATE COMPRESSION IS EASY WITH THE COOLFLEX™

Each Coolflex™ calf garment includes the unique, EZ-Close™ pressure system. The EZ-Close pressure system is positioned laterally, allowing the patient to apply the product without bending forward. The bands of the closure system are made to an optimal length, making them easy to grip and close using minimal hand strength and little dexterity.



The integrated EZ-Close Pressure System has visual indicators for the accurate delivery of 20–30 mmHg, 30–40 mmHg or 40–50 mmHg on each tab



When pairing with a basic liner, a foot piece is required.



When pairing with a transition liner, a foot piece is not required*.

* If a foot piece is not required, please order a transition liner. See page 56.

Standard Calf

20–50 mmHg



Ideal for both venous and lymphatic limb disorders, the garments are available in standard and max sizing and can accommodate limbs that feature disproportionately small ankles with a larger calf.

Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Durable metal clasps, instead of Velcro closure

Includes

- Coolflex™ Standard Calf
- 1 Pair of non-compressive basic Liners per leg

Colour




Black

COOLFLEX™ STANDARD CALF SIZE CHART (A–C1)

(Available in Black Only)

| | Small | Small Max | Medium / Large | Medium / Large Max | X Large / XX Large | |
|-------|----------|---------------|----------------|--------------------|--------------------|---------------|
| C1 | 28–39 cm | 35–50 cm | 34–50 cm | 45–64 cm | 45–64 cm | |
| C | 28–39 cm | 35–50 cm | 34–50 cm | 45–64 cm | 45–64 cm | |
| B | 23–33 cm | 27–38 cm | 29–43 cm | 35–49 cm | 37–54 cm | |
| A | 16–26 cm | 16–26 cm | 23–33 cm | 23–33 cm | 30–42 cm | |
| Left | Regular | 2101-NF-BKR-L | 2101-NF-BMR-L | 2103-NF-BKR-L | 2103-NF-BMR-L | 2105-NF-BKR-L |
| | Tall | 2101-NF-BKT-L | 2101-NF-BMT-L | 2103-NF-BKT-L | 2103-NF-BMT-L | 2105-NF-BKT-L |
| Right | Regular | 2101-NF-BKR-R | 2101-NF-BMR-R | 2103-NF-BKR-R | 2103-NF-BMR-R | 2105-NF-BKR-R |
| | Tall | 2101-NF-BKT-R | 2101-NF-BMT-R | 2103-NF-BKT-R | 2103-NF-BMT-R | 2105-NF-BKT-R |



Standard Foot


Variable compression



Foot compression garments are pictured with the Basic Liner. Basic liner is included with the wrap but not with the footpiece

Features

- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort
- Durable metal clasps, instead of Velcro closure



Foot size measuring Instructions:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure **medial** length from heel to 1st metatarsal head

COOLFLEX™ STANDARD FOOT SIZE CHART (I, J, K)

(Available in Black Only)

| | Small | | Medium / Large | | X Large / XX Large | |
|-------|------------|------------|----------------|------------|--------------------|------------|
| | Regular | Long | Regular | Long | Regular | Long |
| I | 28-39 cm | N/A | 30-43 cm | N/A | 34-49 cm | N/A |
| J | 22-26 cm | N/A | 25-30 cm | N/A | 28-33 cm | N/A |
| K | 14-18 cm | 19-21 cm | 18-20 cm | 20-22 cm | 20-22 cm | 23-25 cm |
| Right | 2101-BTR-R | 2101-BTL-R | 2103-BTR-R | 2103-BTL-R | 2105-BTR-R | 2105-BTL-R |
| Left | 2101-BTR-L | 2101-BTL-L | 2103-BTR-L | 2103-BTL-L | 2105-BTR-L | 2105-BTL-L |

COMPREFLEX™

Wrap Details

ACCURATE, MEASURED THERAPEUTIC COMPRESSION FOR SELF-DIRECTED TREATMENT AND SELF-MANAGED CARE

The Compreflex™ range of products are low profile compression devices designed to deliver moderate to firm compression to the lower extremities. Each Compreflex™ device is made of a soft conforming inelastic material, offering low to minimal stretch for proper fit and the gradient application of compression. The Compreflex™ garments are designed with key features that support the patient, including our patented Accutab® Devices which ensure the accurate application of the prescribed dosage of therapeutic compression, an interior stretch panel for easy donning, a contoured calf for limb shaping, and opposing straps that are rounded for improved fit and patient comfort.

INTELLIGENT DESIGN

Garment features a contoured calf, providing a superior fit

PREMIUM MATERIALS

Soft, conforming material provides a comfortable fit

EASY APPLICATION

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, improving compliance

SUPERIOR QUALITY

Each garment is double stitched to ensure durability



20-50 mmHg



Features & Benefits


- High quality inelastic wrap composed of durable breathable inelastic material with moisture wicking capacities
- Contoured design helps the wrap conform to the natural shape of the leg
- Front stretch panel secures the garment in place, while the Velcro straps assist for easy application
- Ability to customize with different liner and boot options

Conditions

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema

How it works


PROVEN ACCURACY FOR RESTING PRESSURE



Accutab® has been rigorously tested according to established guidelines in the International Compression Club Consensus documents. Testing was conducted using PicoPress, which is recognized for its high accuracy and reproducibility.

Effective working pressure

- Inelastic compression reduces leg edema and improves venous blood flow.



Safe, precise resting pressure

- Under the supervision of a healthcare provider, inelastic compression can be used safely for patients with mild to moderate peripheral arterial disease.
- Accutab® delivers precise dosages for controlled initial resting pressure when donning a Sigvaris Compreflex™ garment.



Compreboot®
Standard

Variable compression
(see page 57 for details)



Basic Liner

Non-Compressive
One pair of Liners
included with each leg
(see page 56 for details)

Standard
CALF & FOOT

20–50 mmHg



For those who require a calf and foot garment

Includes:

- Compreflex™ Wrap
- Non-compressive Basic Liner
- Non attached foot piece

| PRODUCT CODE | | | | | |
|--------------|----------|----------|-----------|----------|----------|
| Size Regular | Black | Beige | Size Tall | Black | Beige |
| Small | 1401-BKR | 1411-BKR | Small | 1401-BKT | 1411-BKT |
| Medium | 1402-BKR | 1412-BKR | Medium | 1402-BKT | 1412-BKT |
| Large | 1403-BKR | 1413-BKR | Large | 1403-BKT | 1413-BKT |
| X Large | 1404-BKR | 1414-BKR | X Large | 1404-BKT | 1414-BKT |
| XX Large | 1405-BKR | 1415-BKR | XX Large | 1405-BKT | 1415-BKT |

Accutab® Devices

(see page 58
for details)



Note: • 1401-BKR, 1411-BKR, 1401-BKT, and 1411-BKT includes a Regular Length Boot

• 1402-BKR, 1412-BKR, 1402-BKT, 1412-BKT, 1403-BKR, 1413-BKR, 1403-BKT, and 1413-BKT includes a Medium/Large Regular Length Boot


• 1404-BKR, 1414-BKR, 1404-BKT, 1414-BKT, 1405-BKR, 1415-BKR, 1405-BKT, and 1415-BKT includes a X-Large/XX-Large Regular Length Boot

Sizing Chart

COMPREFLEX™

CALF SIZING CHART (A–C)

| | Small | Medium | Large | X Large | XX Large |
|---|----------|----------|----------|----------|----------|
| C | 29–39 cm | 34–44 cm | 40–50 cm | 46–56 cm | 54–64 cm |
| A | 19–26 cm | 21–30 cm | 26–36 cm | 31–41 cm | 36–46 cm |



G Length
Regular: 30–33.9 cm
Tall: 34–42 cm

Foot Size:

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure medial length from heel to 1st metatarsal head

Compreflex™ Calf options



Standard
CALF

20-50 mmHg



Transition
CALF

20-50 mmHg



Complete
CALF

20-50 mmHg

For those who require an alternative foot garment

Includes:

- Compreflex™ Wrap
- Non-compressive Basic Liner

Can be matched with alternative foot pieces* (see page 57).

* When paired with a basic liner, a foot piece is required.

For those who require a compressive liner for the foot and ankle

Includes:

- Compreflex™ Wrap
- Transition Liner, which provides 10-15 mmHg in the foot and ankle

Ability to wear regular shoe over liner

For those who require a compressive liner with compression in the foot and ankle, and graduated compression in the calf

Includes:

- Compreflex™ Wrap
- Complete Liner, which provides 10-15 mmHg in the foot and ankle, and graduated compression in the calf

Ability to wear regular shoe over liner

| PRODUCT CODE | | |
|--------------|-------------|-------------|
| Size Regular | Black | Beige |
| Small | 1401-NF-BKR | 1411-NF-BKR |
| Medium | 1402-NF-BKR | 1412-NF-BKR |
| Large | 1403-NF-BKR | 1413-NF-BKR |
| X Large | 1404-NF-BKR | 1414-NF-BKR |
| XX Large | 1405-NF-BKR | 1415-NF-BKR |

| Size Tall | | |
|--------------|-------------|-------------|
| Size Regular | Black | Beige |
| Small | 1401-NF-BKT | 1411-NF-BKT |
| Medium | 1402-NF-BKT | 1412-NF-BKT |
| Large | 1403-NF-BKT | 1413-NF-BKT |
| X Large | 1404-NF-BKT | 1414-NF-BKT |
| XX Large | 1405-NF-BKT | 1415-NF-BKT |

| PRODUCT CODE | | |
|--------------|-------------|-------------|
| Size Regular | Black | Beige |
| Small | 1401-UC-BKR | 1411-UC-BKR |
| Medium | 1402-UC-BKR | 1412-UC-BKR |
| Large | 1403-UC-BKR | 1413-UC-BKR |
| X Large | 1404-UC-BKR | 1414-UC-BKR |
| XX Large | 1405-UC-BKR | 1415-UC-BKR |

| Size Tall | | |
|--------------|-------------|-------------|
| Size Regular | Black | Beige |
| Small | 1401-UC-BKT | 1411-UC-BKT |
| Medium | 1402-UC-BKT | 1412-UC-BKT |
| Large | 1403-UC-BKT | 1413-UC-BKT |
| X Large | 1404-UC-BKT | 1414-UC-BKT |
| XX Large | 1405-UC-BKT | 1415-UC-BKT |

| PRODUCT CODE | | |
|--------------|-------------|-------------|
| Size Regular | Black | Beige |
| Small | 1401-FC-BKR | 1411-FC-BKR |
| Medium | 1402-FC-BKR | 1412-FC-BKR |
| Large | 1403-FC-BKR | 1413-FC-BKR |
| X Large | 1404-FC-BKR | 1414-FC-BKR |
| XX Large | 1405-FC-BKR | 1415-FC-BKR |

| Size Tall | | |
|--------------|-------------|-------------|
| Size Regular | Black | Beige |
| Small | 1401-FC-BKT | 1411-FC-BKT |
| Medium | 1402-FC-BKT | 1412-FC-BKT |
| Large | 1403-FC-BKT | 1413-FC-BKT |
| X Large | 1404-FC-BKT | 1414-FC-BKT |
| XX Large | 1405-FC-BKT | 1415-FC-BKT |

See page 56 for Liner Sizing chart

Liner Options



Basic Liner

No compression

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- drirelease moisture wicking and odour control
- Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Heathered Silver colour only



Transition Liner

10–15 mmHg

- Mild compression in foot and ankle only
 - drirelease moisture wicking and odour control
 - Extra thick sole from toe to heel for added comfort
 - Silver ionized treatment for bacteriostatic properties
 - Available in Black and Beige
- * When ordering a Coolflex standard without a foot piece, please order a transition liner.



Complete Liner

10–15 mmHg

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- drirelease moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties
- Available in Black and Beige

Size: One size fits all

| Pair | 10 Pair Pack |
|------------|--------------|
| 1700 - NK1 | 1700 - NKC |

Included With
Compreflex™ Standard Calf & Foot,
Compreflex™ Standard Calf,
Coolflex™ Standard Calf

PRODUCT CODE (Included with Compreflex™ Transition)

| Size | Pair (Black) | Pair (Beige) |
|----------|--------------|--------------|
| Small | 1601-CK | 1611-CK |
| Medium | 1602-CK | 1612-CK |
| Large | 1603-CK | 1613-CK |
| X Large | 1604-CK | 1614-CK |
| XX Large | 1605-CK | 1615-CK |

PRODUCT CODE (Included with Compreflex™ Complete)

| Size | Pair (Black) | Pair (Beige) |
|----------|--------------|--------------|
| Small | 1601-CK-F1 | 1611-CK-F1 |
| Medium | 1602-CK-F1 | 1612-CK-F1 |
| Large | 1603-CK-F1 | 1613-CK-F1 |
| X Large | 1604-CK-F1 | 1614-CK-F1 |
| XX Large | 1605-CK-F1 | 1615-CK-F1 |

LINERS

TRANSITION AND COMPLETE LINERS SIZING CHART

| Size | Ankle Circ. (A) | Calf Circ. (C) | Women's Shoe | Men's Shoe |
|----------|-----------------|----------------|--------------|------------|
| Small | 19–26 cm | 29–39 cm | <7½ | <6½ |
| Medium | 21–30 cm | 34–44 cm | 7.5–9.5 | 6.5–8.5 |
| Large | 26–36 cm | 40–50 cm | 10–12 | 9–11 |
| X Large | 31–41 cm | 46–56 cm | 12.5–15 | 11.5–14 |
| XX Large | 36–46 cm | 54–64 cm | 12.5–15 | 11.5–14 |

Compreboot® Options



Compreboot®
Lite **FOOT**

Variable
compression

- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild, therapeutic compression



Compreboot®
Standard **FOOT**

Variable
compression

- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape



Compreboot®
Plus **FOOT**

Variable
compression

- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression

| PRODUCT CODE (Available in Black Only) | | |
|--|----------|--|
| Size | Black | |
| Small / Medium / Large | | |
| Regular | 1803-BLR | |
| Long | 1803-BLL | |
| X Large / XX Large | | |
| Regular | 1805-BLR | |
| Long | 1805-BLL | |

| PRODUCT CODE | | |
|---------------------------|----------|----------|
| Size | Black | Beige |
| Small | | |
| Regular | 1801-BTR | 1811-BTR |
| Long | 1801-BTL | 1811-BTL |
| Medium / Large | | |
| Regular | 1803-BTR | 1813-BTR |
| Long | 1803-BTL | 1813-BTL |
| X Large / XX Large | | |
| Regular | 1805-BTR | 1815-BTR |
| Long | 1805-BTL | 1815-BTL |

| PRODUCT CODE (Available in Black Only) | |
|--|----------|
| Size | Black |
| Small | |
| Regular | 1801-BPR |
| Long | 1801-BPL |
| Medium / Large | |
| Regular | 1803-BPR |
| Long | 1803-BPL |
| X Large / XX Large | |
| Regular | 1805-BPR |
| Long | 1805-BPL |

COMPREFLEX® LITE


FOOT SIZING CHART (I, J, K)

| | Small / Medium / Large | | X Large / XX Large | |
|----------|------------------------|----------|--------------------|----------|
| | Regular | Long | Regular | Long |
| I | 28–36 cm | 28–36 cm | ≤44 cm | ≤44 cm |
| J | 22–26 cm | 22–26 cm | ≤33 cm | ≤33 cm |
| K | 14–18 cm | 19–23 cm | 20–22 cm | 25–27 cm |

COMPREFLEX® STANDARD AND PLUS

FOOT SIZING CHART (I, J, K)

| | Small | | Medium / Large | | X Large / XX Large | |
|----------|----------|----------|----------------|----------|--------------------|----------|
| | Regular | Long | Regular | Long | Regular | Long |
| I | 28–36 cm | 28–36 cm | <39 cm | <39 cm | <44 cm | <44 cm |
| J | 22–26 cm | 22–26 cm | <30 cm | <30 cm | <33 cm | <33 cm |
| K | 14–18 cm | 19–23 cm | 18–20 cm | 23–25 cm | 20–22 cm | 25–27 cm |



Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure **medial** length from heel to 1st metatarsal head



Compreflex™ Strap Extender

- Provides up to 10 cm of additional circumferential size range to Compreflex™ calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex™ calf wraps or needing custom size

| | 1 Unit | 10 Unit Pack |
|-------|----------|--------------|
| Black | 1401-SE1 | 1401-SEB |
| Beige | 1411-SE1 | 1411-SEB |



Accutab® Devices

The Accutab® accurately delivers the level of compression prescribed by the health care provider. An easy to use device, the Accutab® ensures the safe delivery of compression in doses of 20–30 mmHg, 30–40 mmHg, or 40–50 mmHg.





Standard Thigh

20–50 mmHg



Features

- The Compreflex™ Standard Thigh is an inelastic compression device for the thigh. The product includes a Belt and Hip Attachment to support the product throughout daily use.
- Pair the Compreflex™ Standard Thigh with the Compreflex™ Standard Knee (sold separately) to create a complete compression system for the upper leg.


Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema

Includes

- Compreflex™ Standard Thigh
- 1 Pair of non-compressive Cotton Liners per leg
- Belt and Hip Attachment
- Accutab® Devices

Colours



COMPREFLEX™ STANDARD THIGH SIZE CHART (D-F) & PRODUCT CODE

| | Small | Medium | Large | X Large |
|----------|----------|----------|----------|----------|
| F | 48–58 cm | 56–66 cm | 64–74 cm | 74–84 cm |
| E | 43–53 cm | 51–61 cm | 58–68 cm | 68–78 cm |
| D | 38–48 cm | 46–56 cm | 53–63 cm | 63–73 cm |

BLACK

| | | | | | |
|-------|---------|------------|------------|------------|------------|
| Left | Regular | 1401-TCR-L | 1402-TCR-L | 1403-TCR-L | 1404-TCR-L |
| | Tall | 1401-TCT-L | 1402-TCT-L | 1403-TCT-L | 1404-TCT-L |
| Right | Regular | 1401-TCR-R | 1402-TCR-R | 1403-TCR-R | 1404-TCR-R |
| | Tall | 1401-TCT-R | 1402-TCT-R | 1403-TCT-R | 1404-TCT-R |

BEIGE

| | | | | | |
|-------|---------|------------|------------|------------|------------|
| Left | Regular | 1411-TCR-L | 1412-TCR-L | 1413-TCR-L | 1414-TCR-L |
| | Tall | 1411-TCT-L | 1412-TCT-L | 1413-TCT-L | 1414-TCT-L |
| Right | Regular | 1411-TCR-R | 1412-TCR-R | 1413-TCR-R | 1414-TCR-R |
| | Tall | 1411-TCT-R | 1412-TCT-R | 1413-TCT-R | 1414-TCT-R |


Standard Knee

20-50 mmHg



Features

- Provides accurate, inelastic compression levels of 20-50 mmHg
- Includes inner donning sleeve for easy application
- Product includes a Butterfly Strap for hard to compress areas



COMPREFLEX™ STANDARD KNEE SIZE CHART (C1-D)

| | Small | Medium | Large | X Large |
|--------------|----------|----------|----------|----------|
| D | 38-48 cm | 46-56 cm | 53-63 cm | 58-73 cm |
| C1 | 29-39 cm | 34-44 cm | 39-49 cm | 44-65 cm |
| Black | 1401-KP | 1402-KP | 1403-KP | 1404-KP |
| Beige | 1411-KP | 1412-KP | 1413-KP | 1414-KP |

*Arrows indicate reference points of D and C1. Length measurement is not required.

COMPREKNEE™

Variable compression



Standard Knee CLOSED PATELLA




Features

- Garment covers the knee, providing full coverage
- Straps provide a full-range of adjustment to create a comfortable fit
- Conforming material for a firm secure fit

Colours



Black Beige



COMPREKNEE™ SIZE CHART (C1-D)

| | Small | Medium | Large / X Large |
|-------|----------|----------|-----------------|
| D | 38–48 cm | 46–56 cm | 53–63 cm |
| C1 | 29–39 cm | 34–44 cm | 39–49 cm |
| Black | 1101-KP | 1102-KP | 1103-KP |
| Beige | 1111-KP | 1112-KP | 1113-KP |

*Arrows indicate reference points of D and C1. Length measurement is not required.

COMPRESHORT®

10-15 mmHg



IDEAL FOR PATIENTS WITH ABDOMINAL OR TRUNCAL EDEMA WHO NEED MILD, SUPPORTIVE COMPRESSION

Product available in short and capri.



Features

The Compresshort® offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema. Available in short and capri lengths, patients with lymphedema and lipedema benefit from the therapeutic compression and support the garment provides.

Ideal for

- Lymphedema
- Lipedema


Colour



Black

| COMPRESHORT® SIZING CHART (Waist, Hip, Midthigh) | | | | PRODUCT CODE | |
|--|------------|------------|------------|---------------------|---------------------|
| | Waist | Hip | Midthigh | Product SKU (Short) | Product SKU (Capri) |
| Small | 69–76 cm | 84–91 cm | 41–46 cm | 1301-CS | 1301-CP |
| Small Max | 76–86 cm | 91–102 cm | 41–46 cm | 1301-CM | 1301-CPM |
| Medium | 76–86 cm | 91–102 cm | 46–53 cm | 1302-CS | 1302-CP |
| Medium Max | 86–97 cm | 102–112 cm | 46–53 cm | 1302-CM | 1302-CPM |
| Large | 86–97 cm | 102–112 cm | 53–61 cm | 1303-CS | 1303-CP |
| Large Max | 97–112 cm | 112–127 cm | 53–61 cm | 1303-CM | 1303-CPM |
| X Large | 97–112 cm | 112–127 cm | 61–71 cm | 1304-CS | 1304-CP |
| X Large Max | 112–132 cm | 127–147 cm | 61–71 cm | 1304-CM | 1304-CPM |
| 2X Large | 112–132 cm | 127–147 cm | 71–81 cm | 1305-CS | 1305-CP |
| 2X Large Max | 132–152 cm | 147–168 cm | 71–81 cm | 1305-CM | 1305-CPM |
| 3X Large | 132–152 cm | 147–168 cm | 81–91 cm | 1306-CS | 1306-CP |
| 3X Large Max | 152–172 cm | 168–188 cm | 81–91 cm | 1306-CM | 1306-CPM |
| 4X Large | 152–172 cm | 168–188 cm | 91–101 cm | 1307-CS | 1307-CP |
| 4X Large Max | 172–193 cm | 188–208 cm | 91–101 cm | 1307-CM | 1307-CPM |
| 5X Large | 172–193 cm | 188–208 cm | 101–112 cm | 1308-CS | 1308-CP |

Note: Available in Black Only



Compression socks and hosiery for everyday wellbeing



Care for Yourself

Why choose **Consumer Health**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

CONSUMER HEALTH PRODUCTS

| Series | Product | Page |
|----------------------|--------------------------|------|
| 146 / 186 | Casual Cotton | 66 |
| 146/T / 186/T | Traveno Travel Socks | 67 |
| 151 / 191 | Sea Island Cotton | 68 |
| 152 / 192 | All-Season Merino Wool | 70 |
| 153 / 193 | Linen | 71 |
| 143 / 183 | Microfiber Shades | 72 |
| 120 | Sheer Fashion | 74 |
| 170 | Soft Silhouette Leggings | 76 |
| 401 | Athletic Recovery Socks | 78 |
| 412V | Performance Sleeves | 79 |
| 421 | Merino Outdoor Socks | 80 |
| 160 | Eversoft Diabetic Socks | 81 |



CASUAL COTTON

146/186

for women & men

Compression Level

15-20 mmHg

Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

Featured colour is Black

SIZING CHART

| | | | | |
|-----|---------------------|---------------------|---------------------|-----------------------|
| 146 | Women's Shoe | 5-7 | 7.5-9.5 | 10-12 |
| | Ankle Circumference | 18-23 cm (7"-9") | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") |
| 186 | Men's Shoe | 6.5-8.5 | 9-11 | 11.5-14 |
| | Ankle Circumference | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") | 21-29 cm (8.5"-11.5") |
| | Size | A | B | C |

STYLES AND COLOURS

Closed Toe = ●



| Calf (C) | White 00 | Navy 10 | Brown 11 | Khaki 30 | Black 99 |
|----------|----------|---------|----------|----------|----------|
| | ● | ● | ● | ● | ● |

Product Code Example:
146CB99

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 146 | C | B | 99 |



Featured colour is Black

TRAVENO TRAVEL SOCKS

146/T - 186/T

for women & men

Compression Level

15-20 mmHg

Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation
- Box cover promotes Traveno as a Travel friendly sock

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

SIZING CHART

| 146/T | Women's Shoe | 5-7 | 7.5-9.5 | 10-12 |
|-------|---------------------|---------------------|---------------------|-----------------------|
| | Ankle Circumference | 18-23 cm (7"-9") | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") |
| 186/T | Men's Shoe | 6.5-8.5 | 9-11 | 11.5-14 |
| | Ankle Circumference | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") | 21-29 cm (8.5"-11.5") |
| Size | A | B | C | |

STYLES AND COLOURS

Closed Toe = ●



| | |
|----------|---|
| Calf (C) | ● |
|----------|---|

Product Code Example: 146CB99/T

| Series | Style | Size | Colour | Travel |
|--------|-------|------|--------|--------|
| 146 | C | B | 99 | /T |

SEA ISLAND COTTON 151/191

for women & men

Compression Level

15-20 mmHg

Also available in 20-30 mmHg, see page 22

Features and Benefits

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free

Featured colour
is Black



15-20 mmHg
Closed Toe Calf Size A in Navy

151CA10



15-20 mmHg
Closed Toe Calf Size B in Brown

191CB11

SIZING CHART

| | | | | |
|-----|---------------------|---------------------|---------------------|-----------------------|
| 151 | Women's Shoe | 5-7 | 7.5-9.5 | 10-12 |
| | Ankle Circumference | 18-23 cm (7"-9") | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") |
| 191 | Men's Shoe | 6.5-8.5 | 9-11 | 11.5-14 |
| | Ankle Circumference | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") | 21-29 cm (8.5"-11.5") |
| | Size | A | B | C |

STYLES AND COLOURS

Closed Toe = ●



| | | | |
|----------|---|---|---|
| Calf (C) | ● | ● | ● |
|----------|---|---|---|

Product Code Example: 151CA99

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 151 | C | A | 99 |



ALL-SEASON MERINO WOOL 152/192

for women & men

Compression Level

15-20 mmHg

Also available in 20-30 mmHg,
see page 24

Features and Benefits

- High quality extra-fine Merino Wool makes this sock soft and breathable
- Moisture wicking capabilities and naturally odour resistant
- Thermoregulating for year-round wear and keeps feet dry
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free

Featured colour
is Black

SIZING CHART

| 152 | Women's Shoe | 5-7 | 7.5-9.5 | 10-12 |
|------|---------------------|---------------------|---------------------|-----------------------|
| | Ankle Circumference | 18-23 cm (7"-9") | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") |
| 192 | Men's Shoe | 6.5-8.5 | 9-11 | 11.5-14 |
| | Ankle Circumference | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") | 21-29 cm (8.5"-11.5") |
| Size | | A | B | C |

STYLES AND COLOURS

Closed Toe = ●



| Calf (C) | NAVY 10 | BROWN 11 | BLACK 99 |
|----------|---------|----------|----------|
| | ● | ● | ● |

Product Code Example: 152CB99

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 152 | C | B | 99 |



Featured colour is Lavender

LINEN 153/193

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 26

Features and Benefits

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free

SIZING CHART

| 153 | Women's Shoe | 5–7 | 7.5–9.5 | 10–12 |
|------|---------------------|---------------------|---------------------|-----------------------|
| | Ankle Circumference | 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 193 | Men's Shoe | 6.5–8.5 | 9–11 | 11.5–14 |
| | Ankle Circumference | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| Size | | A | B | C |

STYLES AND COLOURS

Closed Toe = ●



| Calf (C) | ● | ● | ● | ● |
|----------|---|---|---|---|
| | | | | |

Product Code Example: 153CB207

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 153 | C | B | 207 |

MICROFIBER SHADES 143/183

for women & men

Compression Level

15-20 mmHg

Also available in 20-30 mmHg, see page 30

Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Excellent breathability and moisture management
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
- Latex free

Featured colour is
Royal Blue Argyle



143
Women's
Calf Style

15-20 mmHg
Closed Toe Calf Size B
in Purple Argyle

143CB46



183
Men's
Calf Style

15-20 mmHg
Closed Toe Calf Size C
in Graphite Stripe

183CC17

SIZING CHART

| | | | | |
|------------|---------------------|---------------------|---------------------|-----------------------|
| 143 | Women's Shoe | 5-7 | 7.5-9.5 | 10-12 |
| | Ankle Circumference | 18-23 cm (7"-9") | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") |
| 183 | Men's Shoe | 6.5-8.5 | 9-11 | 11.5-14 |
| | Ankle Circumference | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") | 21-29 cm (8.5"-11.5") |
| | Size | A | B | C |

STYLES AND COLOURS

Women

Closed Toe = ●



| | | | | | | | |
|----------|---|---|---|---|---|---|---|
| Calf (C) | ● | ● | ● | ● | ● | ● | ● |
|----------|---|---|---|---|---|---|---|

Men

Closed Toe = ●



| | | | | | |
|----------|---|---|---|---|---|
| Calf (C) | ● | ● | ● | ● | ● |
|----------|---|---|---|---|---|

Product Code Example: 183CB48

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 183 | C | B | 48 |



Featured colour is Honey

SHEER FASHION 120

for women

Compression Level

15-20 mmHg

Features and Benefits

- Incredible sheer design in a variety of fashionable colours
- Sheer toe and heel area
- Comfortable stretch and soft fabric, excellent for everyday wear

Fiber Content

- 79% Nylon, 21% Spandex
(Content can vary slightly by style)
- Latex free



15-20 mmHg
Closed Toe Calf Size C
in Charcoal

120CC12



15-20 mmHg
Closed Toe Pantyhose Size A
in Black

120PA99

PANTYHOSE SIZING CHART

| cm/ft | Kg/Lbs | 45 / 100 | 48 / 105 | 50 / 110 | 52 / 115 | 54 / 120 | 57 / 125 | 59 / 130 | 61 / 135 | 64 / 140 | 66 / 145 | 68 / 150 | 70 / 155 | 73 / 160 | 75 / 165 | 77 / 170 | 80 / 175 | 82 / 180 | 84 / 185 | 86 / 190 | 89 / 195 | 91 / 200 | 93 / 205 | 95 / 210 | 98 / 215 | 100 / 220 | 102 / 225 | 105 / 230 | | |
|----------------|--------|----------|----------|---------------|---------------------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|-----------|-----------|-----------|--|--|
| 145 cm / 4'9" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 147 cm / 4'10" | | | | SIZE A | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 150 cm / 4'11" | | | | | ANKLE CIRCUMFERENCE | | | | | | | | | | | | | | | | | | | | | | | | | |
| 152 cm / 5'0" | | | | | 19-23 cm (7.5"-9") | | | | | | | | | | | | | | | | | | | | | | | | | |
| 155 cm / 5'1" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 157 cm / 5'2" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 160 cm / 5'3" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 162 cm / 5'4" | | | | | SIZE B | | | | | | | | | | | | | | | | | | | | | | | | | |
| 165 cm / 5'5" | | | | | ANKLE CIRCUMFERENCE | | | | | | | | | | | | | | | | | | | | | | | | | |
| 167 cm / 5'6" | | | | | 19-23 cm (7.5"-9") | | | | | | | | | | | | | | | | | | | | | | | | | |
| 170 cm / 5'7" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 173 cm / 5'8" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 175 cm / 5'9" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 178 cm / 5'10" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 180 cm / 5'11" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 183 cm / 6'0" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

CALF / THIGH SIZING CHART

| 120 | Women's Shoe | 5-7 | 7.5-9.5 | 10-12 |
|-----|---------------------|------------------|---------------------|---------------------|
| | Ankle Circumference | 18-23 cm (7"-9") | 19-25 cm (7.5"-10") | 20-27 cm (8"-10.5") |
| | Size | A | B | C |

Weight Guidelines For Thigh Style: Size A: <135Lbs ; Size B: <180Lbs ; Size C: <230Lbs.

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○



| | CHARCOAL 12 | TAUPE 29 | HONEY 33 | GOLDEN 36 | CAFÉ 73 | BLACK 99 |
|---------------|-------------|----------|----------|-----------|---------|----------|
| Calf (C) | ● | ●○ | ●○ | ●○ | ● | ● |
| Thigh (N) | ● | ● | ● | ● | ● | ● |
| Pantyhose (P) | ● | ● | ● | ● | ● | ● |

Product Code Example: 120NC29

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 120 | N | C | 29 |

Product Code Example: 120PC99

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 120 | P | C | 99 |



SOFT SILHOUETTE LEGGINGS T70

for women

Compression Level

15-20 mmHg

Features and Benefits

- Advanced soft stretch microfiber fabric shapes legs and smooths skin
- Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts

Fiber Content

- 83% Nylon, 17% Spandex
- Latex free



15-20 mmHg
Legging Size B
in Black



Non-binding ankle cuff
stays comfortably in place
for all day wear

170LB99

SIZING CHART


| cm / ft | Kg / Lbs | 45 / 100 | 48 / 105 | 50 / 110 | 52 / 115 | 54 / 120 | 57 / 125 | 59 / 130 | 61 / 135 | 64 / 140 | 66 / 145 | 68 / 150 | 70 / 155 | 73 / 160 | 75 / 165 | 77 / 170 | 80 / 175 | 82 / 180 | 84 / 185 | 86 / 190 | 89 / 195 | 91 / 200 | 93 / 205 | 95 / 210 | 98 / 215 | 100 / 220 | 102 / 225 | 105 / 230 |
|----------------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|-----------|-----------|-----------|
| 147 cm / 4'10" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 150 cm / 4'11" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 152 cm / 5'0" | | | | A | | | | | | | | | | B | | | | | | | | | | | | | | |
| 155 cm / 5'1" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 157 cm / 5'2" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 160 cm / 5'3" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 162 cm / 5'4" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 165 cm / 5'5" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 167 cm / 5'6" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 170 cm / 5'7" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 173 cm / 5'8" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 175 cm / 5'9" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 178 cm / 5'10" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 180 cm / 5'11" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 183 cm / 6'0" | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Guidelines for Ankle Circumference

Size A: 18–23 cm (7"–9"), **Size B & C:** 20–26 cm (8"–10"), **Size D:** 22–29 cm (8.5"–11.5")

STYLES AND COLOURS

Open Foot = ○



Leggings (L)



Product Code Example: 170LB99

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 170 | L | B | 99 |



Featured colour is White



15-20 mmHg
Closed Toe Calf Size M
in Black

401CM99

ATHLETIC RECOVERY SOCKS 401

for women & men

Compression Level

15-20 mmHg

Features and Benefits

- drirelease® yarns provide moisture wicking and odour control
- Comfortable and easy to put on and take off
- Helps reduce lactic acid buildup
- Mild cushioning in sole, heel, and toe area for added comfort
- Smooth, flat toe seam

Fiber Content

- 67% drirelease® (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

SIZING CHART

| | | | | | | |
|-----|---------------------|-----------------------|------------------------|------------------------|--------------------------|--------------------------|
| 401 | Women's Shoe | <7.5 | 7.5-9.5 | 10-12 | 12.5-15 | ≥15.5 |
| | Men's Shoe | <6.5 | 6.5-8.5 | 9-11 | 11.5-14 | 14.5-17 |
| | Ankle Circumference | 18-24 cm (7"-9.5") | 19-26 cm (7.5"-10") | 20-27 cm (8"-10.5") | 21-29 cm (8.5"-11.5") | 21-29 cm (8.5"-11.5") |
| | Size | S | M | L | X | 2 |

STYLES AND COLOURS

Closed Toe = ●

WHITE
00

BLACK
99

| | | |
|----------|---|---|
| Calf (C) | ● | ● |
|----------|---|---|

Product Code Example: 401CM00

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 401 | C | M | 00 |



Featured colour is Purple

PERFORMANCE SLEEVES 412V

for women & men

Compression Level

20-30 mmHg

Features and Benefits

- Stabilizes and supports active muscles
- Ideal for running, triathlons, endurance sports, and gym workouts
- Fabric resists deterioration in salt, chlorine, or fresh water
- Quick drying
- High-visibility reflective transfer

Fiber Content

- 75% Nylon, 25% Spandex
- Latex free

SIZING CHART

| | | | | | |
|------|---------------------|-------------------------|-----------------------------|---------------------------|-------------------------|
| 412V | Ankle Circumference | 18-21.5 cm (7"-8.5") | 21.5-25.5 cm (8.5"-10") | 25.5-29 cm (10"-11.5") | 29-33 cm (11.5"-13") |
| | Calf Circumference | 28-39 cm (11"-15.5") | 31.5-45 cm (12.5"-17.5") | 35.5-51 cm (14"-20") | 40.5-58 cm (16"-23") |
| | Size | S | M | L | X |

STYLES AND COLOURS

| | | | | |
|-----------------|-------------|------------|---------------|-------------|
| Calf Sleeve = ○ | WHITE 00 | BLUE 50 | PURPLE 58* | BLACK 99 |
| Calf Sleeve (V) | ○ | ○ | ○ | ○ |

* Purple is not available in X-Large

Product Code Example: 412VL99

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 412 | V | L | 99 |



Featured colour is Charcoal

MERINO OUTDOOR SOCKS 421

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 38

Features and Benefits

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned, padded soles help prevent blisters
- Mesh flex zones for comfort and ventilation
- Smooth, flat toe seam
- Soft and comfortable
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

SIZING CHART

| | | | | | |
|-----|---------------------|-----------------------|------------------------|------------------------|--------------------------|
| 421 | Women's Shoe | <7.5 | 7.5–9.5 | 10–12 | 12.5–15 |
| | Men's Shoe | <6.5 | 6.5–8.5 | 9–11 | 11.5–14 |
| | Ankle Circumference | 18–24 cm (7"–9.5") | 19–26 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| | Size | S | M | L | X |

STYLES AND COLOURS

Closed Toe = ●

CHARCOAL
12

OLIVE
31

| | | |
|----------|---|---|
| Calf (C) | ● | ● |
|----------|---|---|

Product Code Example: 421CS12

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 421 | C | S | 12 |



EVERSOFT DIABETIC SOCKS 160

for women & men

Compression Level

8–15 mmHg

Features and Benefits

- Padded sole, heel, and toe for comfort and protection
- Smooth, flat toe seam
- Compression-free foot
- drirelease yarns provide moisture wicking and odour control
- Suitable for individuals with limited dexterity or lower tolerance for compression
- Ideal for individuals with sensitive skin

Fiber Content

- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

SIZING CHART

| | | | | | |
|-----|---------------------|-----------------------|------------------------|------------------------|--------------------------|
| 160 | Women's Shoe | <7.5 | 7.5–9.5 | 10–12 | 12.5–15 |
| | Men's Shoe | <6.5 | 6.5–8.5 | 9–11 | 11.5–14 |
| | Ankle Circumference | 18–24 cm (7"–9.5") | 19–26 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| | Size | S | M | L | X |

STYLES AND COLOURS

Closed Toe = ●

WHITE
00

| | |
|----------|---|
| Calf (C) | ● |
|----------|---|

Product Code Example: 160CS00

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 160 | C | S | 00 |

Donning Aids and Product Care Accessories



Make Life Easier

Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Extra liners for added comfort

ACCESSORIES

| Product | Page |
|------------------------------|------|
| Gloves | 84 |
| Cotton Underliner | 85 |
| Washing Solution | 85 |
| It Stays Roll-On Adhesive | 85 |
| Sim-Slide | 86 |
| Magnide On/Off | 86 |
| Sigvaris Doff N' Donner | 87 |
| Sigvaris Doff N' Donner Cone | 87 |
| Melany | 87 |

GLOVES



Rubber Gloves*

Ridged - The practical and protective way to apply and remove your compression stockings.



Latex-Free Gloves

Gripping donning gloves recommended for use in a latex-free environment.



Textile Gloves*

These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.

PACKAGE: 12 PER PACK

| SIZE | PRODUCT CODE |
|--------------------|--------------|
| Small | 591R400S |
| Medium | 591R400M |
| Large | 591R400L |
| X-Large | 591R400X |
| ASST'D 4X S/M/L | 591R400A |

PACKAGE: 12 PER PACK

| SIZE | PRODUCT CODE |
|---------|--------------|
| Small | 592R400S |
| Medium | 592R400M |
| Large | 592R400L |
| X-Large | 592R400X |

PACKAGE: 12 PER PACK

| SIZE | PRODUCT CODE |
|---------------------|--------------|
| X-Small | 593R400XS |
| Small | 593R400S |
| Medium | 593R400M |
| Large | 593R400L |
| X-Large | 593R400X |
| ASST'D 4x XS/S/M | 593R400A |

RUBBER GLOVE AND LATEX-FREE GLOVE SIZING CHART

Measure Circumference of the Hand at the Palm

| Hand Circumference | 16.5–18 cm (6.5"–7") | 19–20.5 cm (7.5"–8") | 21.5–23 cm (8.5"–9") | 24–25.5 cm (9.5"–10") |
|--------------------|----------------------|----------------------|----------------------|-----------------------|
| | Small | Medium | Large | X-Large |

TEXTILE GLOVE SIZING CHART

Measure Length and Width of the Hand

| Hand Length** | 20.5–21 cm | 21.5–22 cm | 22.5–23 cm | 23.5–24 cm | 24.5–25 cm |
|---------------|------------|------------|------------|------------|------------|
| Hand Width | 7–7.5 cm | 7.5–8 cm | 8–8.5 cm | 9–9.5 cm | 10–10.5 cm |
| | X-Small | Small | Medium | Large | X-Large |

* **Caution:** Rubber Gloves and Textile Gloves contain natural rubber latex that may cause allergic reaction.

** **Hand Length:** measured from the tip of the longest finger to the crease under the palm.

LINERS



Cotton Underliner

10–15 mmHg

Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

For sizing, please refer to the Natural Rubber Calf sizing on page 43.

PACKAGE: 2 PER PACK

| SIZE | PRODUCT CODE | NATURAL RUBBER SIZING |
|------------------|--------------|-----------------------|
| Small / Short | 588USS | S1/S3 |
| Small / Long | 588USL | S2/S4 |
| Medium / Short | 588UMS | M1/M3 |
| Medium / Long | 588UML | M2/M4 |
| Large-XL / Short | 588ULS | L1/L3 - X1/X3 |
| Large-XL / Long | 588ULL | L2/L4 - X2/X4 |

WASHING SOLUTION

Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!



PACKAGE OPTIONS

| CONTAINS | PRODUCT CODE |
|-----------------------------|--------------|
| Six 16 oz. Bottles (473ml) | 586W900 |
| One 16 oz. Bottle (473ml) | 586W950 |
| Twelve 2 oz. Bottles (60ml) | 586W975 |

IT STAYS!

Roll-on body adhesive

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin.

IT STAYS! is pliable and moves with the body.



| SIZE | PRODUCT CODE |
|-----------------------------|--------------|
| One 2 oz Bottle (60ml) | 5811350 |
| Twelve 2 oz. Bottles (60ml) | 5811300 |

DONNING AIDS



Sim-Slide

(For open toe Styles)

| SIZE | PRODUCT CODE |
|---------|--------------|
| Small | 587SLIDESM |
| Medium | 587SLIDEMD |
| Large | 587SLIDELG |
| X-Large | 587SLIDEXL |

SIZING CHART For Sim-Slide

| | | | | |
|-------------------|----------|----------|----------|-----------|
| Women's Shoe Size | ≤4-6 | 6.5-8 | 8.5-9.5 | ≥10 |
| Men's Shoe Size | ≤3-4.5 | 5-6.5 | 7-8 | ≥8.5 |
| Size | S | M | L | XL |



Magnide On/Off

(For closed toe Styles)

| SIZE | PRODUCT CODE |
|---------|--------------|
| Medium | 587MAGMD |
| Large | 587MAGLG |
| X-Large | 587MAGXL |

SIZING CHART For Magnide On/Off

| | | | |
|-------------------|----------|----------|-----------|
| Women's Shoe Size | ≤8.5 | 9-12 | ≥ 12.5 |
| Men's Shoe Size | ≤7.5 | 8-11 | ≥ 11.5 |
| Size | M | L | XL |

* XL suitable for extra-wide feet

DONNING AIDS



SIGVARIS **Doff N' Donner**

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.

| SIZE | PRODUCT CODE |
|----------|--------------|
| One Size | 583D100 |



SIGVARIS **Doff N' Donner Cone**

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

| SIZE | PRODUCT CODE |
|----------|--------------|
| One Size | 583D200 |



MELANY

MELANY from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.*

| SIZE | CALF CIRCUMFERENCE | PRODUCT CODE |
|---------|--------------------|--------------|
| X-Small | 25–30 cm | 579SXS |
| Small | 30–38 cm | 579SSM |
| Medium | 38–46 cm | 579SMD |
| Large | 46–56 cm | 579SLG |

* Individual components not sold separately

Learn more



Help yourself

How can **resources** help?

- Learn about the levels of compression for various indications to better help your clients
- Understand the terms and icons that we use
- Follow specific instructions on how to measure for the perfect fit

RESOURCES

| Section | Page |
|------------------------|------|
| Symptoms & Indications | 90 |
| Glossary of Icons | 91 |
| Glossary of Terms | 92 |
| How to Measure | 93 |
| How to Order | 94 |

| SYMPTOMS / INDICATIONS | COMPRESSION LEVEL AT THE ANKLE |
|---|---------------------------------|
| <ul style="list-style-type: none"> □ Tired, achy legs □ Occupational or evening edema □ Leg discomfort from long hours of standing or sitting □ Prevent varicosities & swelling during pregnancy □ Preventative for legs predisposed to risk □ Reduce swelling during travel □ Prevent DVT during travel | <p>15–20 mmHg</p> |
| <ul style="list-style-type: none"> □ Heavy, fatigued, aching legs □ Mild edema in lower extremities □ Mild varicosities with minimal edema □ Mild varicosities during pregnancy □ Lymphedema upper extremities □ Post-procedure of small veins □ Prevention or management of DVT <p><i>Note: A prescription is recommended for stockings with the indications of 20 mmHg or higher</i></p> | <p>20–30 mmHg</p> |
| <ul style="list-style-type: none"> □ Moderate varicosities with mild to moderate edema (also during pregnancy) □ Post-procedure of larger veins to maintain treatment success □ CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations) □ Management and treatment of DVT or Post-Thrombotic Syndrome □ Superficial phlebitis □ Lymphedema after decongestant therapy to maintain reduction □ Orthostatic / postural hypotension □ Post-fracture, post-traumatic edema | <p>30–40 mmHg</p> |
| <ul style="list-style-type: none"> □ Lymphedema after decongestant therapy to maintain reduction □ Severe Post-Thrombotic Syndrome □ Severe varicosities and/or edema | <p>40 mmHg or higher</p> |
| <ul style="list-style-type: none"> □ Severe arterial insufficiency, intermittent claudication, ischemia □ Uncontrolled congestive heart failure | <p>Contraindications</p> |
| <ul style="list-style-type: none"> □ Skin sensitivities or allergies □ Advanced neuropathy with sensory loss in leg or foot □ Diabetes with arterial compromise □ Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician □ Eczema, hypodermatitis, dermatosclerosis, dermatitis | <p>Cautions</p> |
| <p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p> | |
| <p><i>No liability accepted for non-observance of contraindications and cautions.</i> Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p> | |



Sensitive Skin



Thermoregulating



Breathable



Moisture Wicking



Soft



Durable



Cotton



Wool



Microfiber



Natural Rubber



Largest Size Range



Sheer



Medium Sheer



Opaque



Stylish



Stylish



Elegant



Colours



Patterns



Sensinno® Grip-top Band



Anatomically Contoured



Cushioned Sole



Anti-bacterial



Hypoallergenic Natural fibers



Moisture Management

| TERM | DEFINITION | SIGVARIS QUALITY STANDARD |
|---|--|--|
| CHRONIC VENOUS DISEASE (CVD) | A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed). | SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD. |
| DOUBLE-COVERED YARN | Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort. | Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort. |
| DEEP VEIN THROMBOSIS (DVT) | A blood clot in a deep vein, most commonly in the lower leg or thigh. | SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening. |
| GRADUATED COMPRESSION VS. ANTI-EMBOLISM | Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-ambulatory patients. | SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation. |
| HYPOALLERGENIC | Designed to minimize the chance of an allergic reaction. | SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin. |
| ISO CERTIFIED | SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place. | SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America. |
| MICROFIBER | An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural fiber fabric. This fabric has excellent breathability. | Used to create products with a soft and comfortable feel. Find this fiber in Soft Opaque 840, Microfiber Patterns 830, Microfiber Shades 143/183, and Microfiber 820. |
| mmHg (MILLIMETERS OF MERCURY) | Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure. | SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement. |
| MERINO WOOL | Merino is the softest and finest wool available in the world and is inherently itch-free. | SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in Thermoregulating Wool 422, Merino Outdoor Socks 421, Merino Wool 240, and All-Season Merino Wool 152/192. |
| SEA ISLAND COTTON | Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness. | SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our Sea Island 220 and Sea Island Cotton 151/191. |
| SENSINNOV® | A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day. | This top band is featured in Sheer 780, Soft Opaque 840, and Microfiber 820 to help keep thigh-high stockings in place all day long. |
| SUPIMA COTTON | Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber. | SIGVARIS GROUP uses this premium cotton yarn in the Cotton 230, Sea Island 220, Sea Island Cotton 151/191, Casual Cotton 146/186, and Cushioned Cotton 360, making them comfortable and ideal for those with sensitive skin. |

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE

Open/Closed Toe

| | |
|----------|---|
| 1 | Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement. |
| 2 | Measure the calf circumference at the fullest part of the calf. |
| 3 | Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg. |



1 Ankle circumference directly above ankle bone



2 Calf circumference at fullest part of the calf


PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe

| | |
|----------|--|
| 1 | Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be your first measurement. |
| 2 | Measure the calf circumference at the fullest part of the calf. |
| 3 | Skip Step 3 (but refer to General Guidelines #6 below). |
| 4 | Measure the thigh at its widest circumference. |
| 5 | Determine the leg length by measuring from the floor up to the gluteal fold. |
| 6 | Determine the hip circumference by measuring the widest part of the hips. |



3 Calf length from the floor to the right-angle bend of the knee



4 Widest circumference of the thigh

GENERAL GUIDELINES

| | |
|----------|--|
| 1 | Measurements should always be taken as early in the morning as possible, when swelling is minimal. |
| 2 | Measure against bare skin, not over clothing (Steps 1-5). |
| 3 | The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples). |
| 4 | Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting. |
| 5 | Length measurements should be taken from the ground up. |
| 6 | Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles. |
| 7 | Record shoe size; it may be helpful in determining the best possible fit. |



5 Leg length from the floor up to the gluteal fold



6 Circumference at the hip

HOW TO ORDER MEDICAL AND SPECIALTY PRODUCTS



Product Code Example: 232CMLM99 (page 8)

| Series Number | Compression Level | Product Style | Size | Gender | Colour |
|--|--|---|--|--------------------------------------|---|
| 23 = Cotton <i>(See individual pages, use first 2 digits of series number)</i> | 1 = 15–20 mmHg 2 = 20–30 mmHg 3 = 30–40 mmHg 4 = 40–50 mmHg* 5 = 50–60 mmHg* | C = Knee-high or Calf <i>(See individual pages for available styles)</i> | ML = Medium Long <i>(See individual pages for sizing)</i> | M = Men W = Women O = Open toe | 99 = Black <i>(See individual pages for available colours)</i> |
| 23 | 2 | C | ML | M | 99 |
| * 40–50 mmHg and 50–60 mmHg available in Natural Rubber only | | | | | |
| To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only): add “/P” after the colour in the ordering code | | | | | |
| To order Calf Style with Grip-top: add “/S” after the colour in the ordering code <i>(see product pages for availability)</i> | | | | | |
| To order Thigh with Waist Attachment: add “/L” to indicate the Left leg or “/R” to indicate the Right leg after the colour in the ordering code <i>(860W and 500W only)</i> | | | | | |
| To order Motion High Tech and Motion Thermoregulating Wool: The letter for Gender is not used in the ordering code. Example: 412CSM59 = Motion High Tech, 20–30 mmHg, Closed toe Calf, Size Small Medium, Colour Red | | | | | |

HOW TO ORDER CONSUMER HEALTH PRODUCTS



Product Code Examples:

151CA99 (page 69)

191CB99 (page 69)

| Series Number | Product Style | Size | Colour |
|--|--|--|---|
| 151/191 = SEA ISLAND COTTON First 3 digits indicate gender 151 = woman 191 = man <i>(See individual pages for series number)</i> | C = Calf <i>(See individual pages for available styles)</i> | According to shoe size A= women's shoe sizes 5-7 B = men's shoe sizes 9-11 <i>(See individual pages for sizing)</i> | 99 = Black <i>(See individual pages for available colours)</i> |
| 151 | C | A | 99 |
| 191 | C | B | 99 |
| <p>ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFIT DIABETIC SOCKS are UNISEX: the first 3 digits are the same for both genders.</p> <p>Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20 mmHg, Closed toe Calf, Size Large, Colour Olive</p> | | | |

Sigvaris, SIGVARIS GROUP, PFS², Sensinno, Compreflex, Compreboot, and Accutab are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2021 SIGVARIS, Inc. Coolflex, Spacer Fabric TR, and CoverOne are trademarks of SIGVARIS AG, CH-9014, St. Gallen/Switzerland, in many countries worldwide. Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. NILIT is a registered trademark of Nilit. All rights reserved. *drire/lease* is copyright ©2021 *drire/lease*. *drire/lease* is an Optimer Brands company. All rights reserved. Sim-Slide and Magnide are registered trademarks of Arion Holding B.V.

mobilis

by SIGVARIS GROUP

**SOFT ORTHOPEDIC
SUPPORTS AND
BRACES**



INDEX

Intro

98

Mode of action 99
How to measure 99

Elbow

100

EPICARE - Elbow Support 101
EPIBAND - Elbow Strap 101
EPIACTIVE - Elbow Brace 102

Wrist and thumb

103

MANUCARE - Wrist Support 104
MANUWRAP - Wrist and Thumb Brace 104
MANUACTIVE - Wrist Brace 105
MANUSUPPORT - Wrist Brace 105
MANUPLUS SUPPORT - Wrist and Thumb Stabilizer 106
RHIZOSUPPORT - Thumb Brace 106

Back

107

LUMBOCARE - Back Support 108
LUMBOSUPPORT - Lumbar Brace 108
LUMBOSUPPORT EVO - Lumbar Brace 109
LUMBOINTENSE - Lumbar Brace 109

Knee

110

GENUCARE - Knee Support 111
GENUWRAP - Knee Wrap 111
GENUBAND - Knee Strap 112
GENUACTIVE STAB - Stabilizing Knee Brace 112
GENUACTIVE PAD - Patellar Knee Brace 113

Ankle

114

MALLEOCARE - Ankle Support 115
MALLEOACTIVE PAD - Ankle Brace 115
MALLEOACTIVE BAND - Ankle Brace 116
MALLEOSUPPORT - Ankle Brace 116
MALLEOINTENSE - Ankle Stabilizer 117

How to order Mobilis

117

MOBILIS: SOFT ORTHOPEDIC SUPPORTS AND BRACES

A collection of high-quality supports and braces from SIGVARIS GROUP designed to provide support, stability, and pain relief while allowing individuals to stay active.




Mobilis product names are made of two parts:

1 Body part where the product is used


- Elbow → **EPI**
- Hand → **MANU**
- Thumb → **RHIZO**
- Back → **LUMBO**
- Knee → **GENU**
- Ankle → **MALLEO**

2 Support level:




Example **1** + **2** = **MALLEOSUPPORT**

MODE OF ACTION




Compression effect

- Decreases swelling and pressure on joints
- Provides pain relief
- Improves mobility



Massage effect with padding

- Targeted swelling reduction
- Improves blood circulation
- Aids healing



Joint stabilization

- Rigid elements stabilize the joint and improve muscle activity



Proprioception

- Aids in perception of movements and body position




Warming effect


- Increases comfort and may provide pain relief

HOW TO MEASURE

Elbow



Measure the arm circumference at 3 cm below the elbow (two fingers)




Measure the elbow circumference with 30 degrees of inclination

Wrist and thumb




Measure wrist circumference

Back



Measure waist circumference at belly button

Knee



Measure knee circumference with 30 degrees of inclination



Measure thigh circumference at 15 cm above the knee

Ankle



Measure ankle circumference above ankle bone



Elbow



| | EPI CARE | EPI BAND | EPI ACTIVE |
|--|----------------------------------|------------------------------------|------------------------------------|
| INDICATIONS | | | |
| Arthritis | ○ | | ○ |
| Joint fluid and swelling | ○ | | ○ |
| Chronic irritation/ Tendomyopathy | ○ | | ○ |
| Post trauma, Post surgery | ○ | | ○ |
| Tennis/golfer's elbow (Epicondylopathy) | ○ | ○ | ○ |
| SUPPORT LEVEL | Support level mild | Support level medium | Support level medium |

EPI CARE

Elbow Support



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)



Contains natural rubber latex that may cause allergic reaction

ACTION MODE




FEATURES

- Mild support
- Maximum comfort
- Discreet look
- Can be worn on either elbow
- Optimal adhesion due to silicone band

COMPOSITION

64% Polyamide, 27% Elastodiene, 5% Elastane, 4% Silicone



Measure the elbow circumference with 30 degrees of inclination



| Size | Circumference | SKU |
|----------------|---------------|------------|
| Size XS | 20-23 cm | MECE-XS-66 |
| Size S | 23-26 cm | MECE-SM-66 |
| Size M | 26-29 cm | MECE-MD-66 |
| Size L | 29-32 cm | MECE-LG-66 |
| Size XL | 32-35 cm | MECE-XL-66 |

EPI BAND

Elbow Strap



INDICATIONS

Tennis/golfer's elbow (epicondylopathy)

ACTION MODE




FEATURES

- Adjustable strap to modify support level
- Silicone pad provides relief for tendon
- Can be worn on either elbow

COMPOSITION

89% Polyamide, 11% Elastane



Measure the arm circumference at 3 cm below the elbow (two fingers)



| Size | Circumference | SKU |
|----------|---------------|---------|
| One size | 20 - 35 cm | MEBE-91 |

EPIACTIVE

Elbow Brace



Contains natural rubber latex that may cause allergic reaction



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

ACTION MODE



Compression



Massaging effect




Proprioception

FEATURES

- Elbow brace with silicone inserts
- Removable elastic strap to adjust support
- Anatomical shape and 3D knitting for maximum comfort
- Can be worn on either elbow

COMPOSITION

56% Thermoplastic elastomer, 32% Polyamide, 7% Elastodiene, 3% Elastane, 2% Polyester



Measure the elbow circumference with 30 degrees of inclination

| Size | Circumference | SKU |
|----------------|---------------|------------|
| Size XS | 20 - 23 cm | MEAE-XS-91 |
| Size S | 23 - 26 cm | MEAE-SM-91 |
| Size M | 26 - 29 cm | MEAE-MD-91 |
| Size L | 29 - 32 cm | MEAE-LG-91 |
| Size XL | 32 - 35 cm | MEAE-XL-91 |



Click here to watch a video



Wrist & Thumb



| | MANU CARE | MANU WRAP | MANU ACTIVE | MANU SUPPORT | MANU PLUS SUPPORT | RHIZO SUPPORT |
|---|--------------------|------------------------|----------------------|----------------------|----------------------|----------------------|
| INDICATIONS | | | | | | |
| Arthritis | ○ | ○ | ○ | ○ | ○ | ○ |
| Chronic irritation/ Tendomyopathy | ○ | ○ | ○ | ○ | ○ | ○ |
| Carpal tunnel syndrome | | | ○ | ○ | ○ | |
| Instability | | | ○ | ○ | ○ | ○ |
| Post trauma, Post surgery | | ○ | ○ | ○ | ○ | |
| Rhizarthrosis | | | | | ○ | ○ |
| Skier's thumb/Torn thumb ligament (Lateral metacarpophalangeal ligament lesion) | | | | | ○ | ○ |
| SUPPORT LEVEL | Support level mild | Support level moderate | Support level medium | Support level strong | Support level strong | Support level strong |

MANUCARE

Wrist Support



INDICATIONS

Arthritis, chronic irritation/tendomyopathy

ACTION MODE



Compression



Warming effect



Proprioception

FEATURES

- Mild support
- Soft compression
- Ultimate comfort
- Discreet look
- Can be worn on either wrist

COMPOSITION

66% Polyamide, 29% Elastodiene, 5% Elastane



Contains natural rubber latex that may cause allergic reaction



Measure wrist circumference

| Size | Circumference | SKU |
|----------------|---------------|------------|
| Size XS | 13 - 15 cm | MMCW-XS-66 |
| Size S | 15 - 17 cm | MMCW-SM-66 |
| Size M | 17 - 19 cm | MMCW-MD-66 |
| Size L | 19 - 21 cm | MMCW-LG-66 |
| Size XL | 21 - 23 cm | MMCW-XL-66 |

MANUWRAP

Wrist and Thumb Brace



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery

ACTION MODE



Movement control



Stabilization

FEATURES

- Moderate adjustable support level
- Breathable fabric
- Can be worn on either hand

COMPOSITION

96% Polyamide, 4% Elastane



Measure wrist circumference



| Size | Circumference | SKU |
|-------------------|---------------|---------------|
| Sizes S/M | 13 - 17 cm | MMWW-SM/MD-91 |
| Sizes L/XL | 17 - 23 cm | MMWW-LG/XL-91 |


MANUACTIVE

Wrist Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

ACTION MODE



FEATURES

- Brace with silicone inserts around cubitus bone
- Metal rod to limit wrist movement
- Elastic strap allows support level adjustment
- Anatomically shaped 3D knit for comfort and fit
- Specific for left or right hand

COMPOSITION

56% Polyamide, 25% Thermoplastic Elastomer, 9% Elastodiene, 5% Polyester, 5% Elastane



Measure wrist circumference



Support level
medium



Contains natural rubber latex that may cause allergic reaction



[Click here to watch a video](#)



| Size | Circumference | SKU Left Hand | SKU Right Hand |
|----------------|---------------|---------------|----------------|
| Size XS | 13 - 15 cm | MMAW-XS-L-91 | MMAW-XS-R-91 |
| Size S | 15 - 17 cm | MMAW-SM-L-91 | MMAW-SM-R-91 |
| Size M | 17 - 19 cm | MMAW-MD-L-91 | MMAW-MD-R-91 |
| Size L | 19 - 21 cm | MMAW-LG-L-91 | MMAW-LG-R-91 |
| Size XL | 21 - 23 cm | MMAW-XL-L-91 | MMAW-XL-R-91 |


MANUSUPPORT

Wrist Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

ACTION MODE



FEATURES

- Immobilization/stabilization brace: rigid palm structure, external and internal rods to restrain wrist movements
- Three hook-and-loop straps to adjust the fit
- Can be worn on either wrist

COMPOSITION

92% Polyamide, 5% Elastane, 3% Cotton



Support level
strong



[Click here to watch a video](#)



Measure wrist circumference

| Size | Circumference | SKU |
|-------------------|---------------|----------------|
| Sizes S/M | 13 - 17 cm | MMSWS-SM/MD-91 |
| Sizes L/XL | 17 - 23 cm | MMSWS-LG/XL-91 |

MANUPLUS SUPPORT


Wrist & Thumb Stabilizer



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

ACTION MODE



FEATURES

- Immobilization/stabilization brace for wrist and thumb
- The two hook-and-loop straps can be clamped by laces to adjust the fit
- Specific for right or left hand

COMPOSITION

50% Polyester, 46% Polyamide, 4% Polyurethane



| Size | Circumference | SKU Left Hand | SKU Right Hand |
|----------------|---------------|---------------|----------------|
| Size XS | 13 - 15 cm | MMPSW-XS-L-91 | MMPSW-XS-R-91 |
| Size S | 15 - 17 cm | MMPSW-SM-L-91 | MMPSW-SM-R-91 |
| Size M | 17 - 19 cm | MMPSW-MD-L-91 | MMPSW-MD-R-91 |
| Size L | 19 - 21 cm | MMPSW-LG-L-91 | MMPSW-LG-R-91 |
| Size XL | 21 - 23 cm | MMPSW-XL-L-91 | MMPSW-XL-R-91 |

RHIZOSUPPORT


Thumb Brace



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, instability, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

ACTION MODE



FEATURES


- Immobilization/stabilization brace for thumb due to external and internal rods
- Two hook-and-loop straps to adjust the fit
- Can be worn on either hand

COMPOSITION

92% Polyamide, 8% Elastane




Click here to watch a video



| Size | Circumference | SKU |
|-------------------|---------------|---------------|
| Sizes S/M | 13 - 17 cm | MRST-SM/MD-91 |
| Sizes L/XL | 17 - 23 cm | MRST-LG/XL-91 |



Back



LUMBO CARE **LUMBO SUPPORT** **LUMBO SUPPORT EVO** **LUMBO INTENSE**

INDICATIONS

| | | | | |
|---|-------------------------------|---------------------------------|---------------------------------|---------------------------------------|
| Muscular tension | ○ | ○ | ○ | ○ |
| Lower back pain (Lumbar syndrome) | ○ | ○ | ○ | ○ |
| Muscular weakness and loss (degeneration) | | ○ | ○ | ○ |
| Post surgery | | ○ | ○ | ○ |
| Sciatica | | ○ | ○ | ○ |
| SUPPORT LEVEL | Support level mild | Support level strong | Support level strong | Support level extra strong |

LUMBOCARE


Back Support



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome)

ACTION MODE




FEATURES

- Adjustable belt provides mild support
- Comfortable material that provides a warming effect
- Can be worn against the skin

COMPOSITION

65% Cotton, 31% Polyester, 4% Elastodiene



| Size | Circumference | SKU |
|---------------|---------------|------------|
| Size S | 70 - 95 cm | MLCB-SM-66 |
| Size M | 95 - 115 cm | MLCB-MD-66 |
| Size L | 115 - 135 cm | MLCB-LG-66 |

LUMBOSUPPORT


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



FEATURES


- Four stabilizing rods to support and secure
- Adjustable straps for a precise fit
- To be worn over clothing or underwear
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

54% Polyamide, 21% Elastane, 20% Polypropylene, 5% Polyester



Click here to watch a video



| Size | Circumference | SKU Short (S) | SKU Tall (T) |
|---------------|---------------|---------------|--------------|
| Size S | 70 - 95 cm | MLSB-SM-S-91 | MLSB-SM-T-91 |
| Size M | 95 - 115 cm | MLSB-MD-S-91 | MLSB-MD-T-91 |
| Size L | 115 - 135 cm | MLSB-LG-S-91 | MLSB-LG-T-91 |

LUMBOSUPPORT EVO


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



Click here to watch a video

FEATURES


- Lumbar spine stabilizing belt with four posterior rods
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

44% Polyamide, 23% Polyester, 8% Viscose, 7% Elastodiene, 6% Ethylene vinyl acetate, 5% Polyvinyl chloride, 4% Polyurethane, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



| Size | Circumference | SKU Short (S) | SKU Tall (T) |
|-----------------|---------------|---------------|---------------|
| Size S | 56 - 71 cm | MLSEB-SM-S-99 | MLSEB-SM-T-99 |
| Size M | 68 - 84 cm | MLSEB-MD-S-99 | MLSEB-MD-T-99 |
| Size L | 81 - 96 cm | MLSEB-LG-S-99 | MLSEB-LG-T-99 |
| Size XL | 94 - 109 cm | MLSEB-XL-S-99 | MLSEB-XL-T-99 |
| Size XXL | 107 - 122 cm | MLSEB-2X-S-99 | MLSEB-2X-T-99 |

LUMBOINTENSE


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



Click here to watch a video

FEATURES


- Stabilizing belt with four posterior rods provides stronger pressure on the lumbar spine
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

47% Polyamide, 27% Polyester, 2% Viscose, 4% Elastodiene, 6% Ethylene-acetate, 5% Polyvinyl chloride, 4% Polyurethane, 2% Polypropylene, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



| Size | Circumference | SKU Short (S) | SKU Tall (T) |
|-----------------|---------------|---------------|--------------|
| Size S | 56 - 71 cm | MLIB-SM-S-99 | MLIB-SM-T-99 |
| Size M | 68 - 84 cm | MLIB-MD-S-99 | MLIB-MD-T-99 |
| Size L | 81 - 96 cm | MLIB-LG-S-99 | MLIB-LG-T-99 |
| Size XL | 94 - 109 cm | MLIB-XL-S-99 | MLIB-XL-T-99 |
| Size XXL | 107 - 122 cm | MLIB-2X-S-99 | MLIB-2X-T-99 |

Knee



| | GENU CARE | GENU WRAP | GENU BAND | GENU ACTIVE STAB | GENU ACTIVE PAD |
|--|--------------------|------------------------|----------------------|----------------------|----------------------|
| INDICATIONS | | | | | |
| Arthritis | ○ | ○ | | ○ | ○ |
| Joint fluid and swelling | ○ | | | ○ | ○ |
| Chronic irritation/ Tendomyopathy | ○ | ○ | ○ | ○ | ○ |
| Post trauma, Post surgery | ○ | ○ | | ○ | ○ |
| Kneecap (patellar) pain | | ○ | ○ | ○ | ○ |
| Deep knee pain (Femoro-patellar syndrome) | | ○ | ○ | ○ | ○ |
| Instability, Hyperlaxity, Meniscus lesion, Gonarthrosis | | | | ○ | |
| SUPPORT LEVEL | Support level mild | Support level moderate | Support level medium | Support level medium | Support level medium |

GENUCARE

Knee Support



Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery

ACTION MODE




FEATURES

- Mild support
- Optimal adhesion due to silicone band
- Discreet look
- Can be worn on either knee

COMPOSITION

64% Polyamide, 28% Elastodiene
5% Elastane, 3% Silicone



Measure circumference below the kneecap with 30 degrees of inclination

| Size | Circumference | SKU |
|----------------|---------------|------------|
| Size XS | 30 - 34 cm | MGCK-XS-66 |
| Size S | 34 - 38 cm | MGCK-SM-66 |
| Size M | 38 - 42 cm | MGCK-MD-66 |
| Size L | 42 - 46 cm | MGCK-LG-66 |
| Size XL | 46 - 50 cm | MGCK-XL-66 |



GENUWRAP

Knee Wrap



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Moderate support knee wrap with three straps for optimum fit
- Adjustable size
- Breathe-O-Prene material allows skin to breathe
- Can be worn on either knee

COMPOSITION

94% Polyamide, 6% Elastane



Measure circumference below the kneecap with 30 degrees of inclination

| Size | Circumference | SKU |
|------------------|---------------|---------------|
| Size S/M | 30 - 42 cm | MGWK-SM/MD-91 |
| Size L/XL | 42 - 55 cm | MGWK-LG/XL-91 |



GENUBAND

Knee Strap



INDICATIONS

chronic irritation/tendomyopathy, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Adjustable knee strap to modify support level
- Silicone insert is located on kneecap tendon
- Can be worn on either knee
- Anatomical shape for extra comfort

COMPOSITION

89% Polyamide, 11% Elastane



| Size | Circumference | SKU |
|----------|---------------|---------|
| One Size | 30 - 55 cm | MGBK-91 |

GENUACTIVE STAB

Stabilizing Knee Brace



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome), instability, hyperlaxity, meniscus lesion, gonarthrosis

ACTION MODE



FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap with three Velcro straps
- The anatomically shaped fabric is compressive
- Wearable on both knees

COMPOSITION


8% Polyamide, 23% Elastomer, 15% Elastodiene, 6% Polyester, 5% Elastane, 2% Polyvinyl chloride, 1% Silicone



Contains natural rubber latex that may cause allergic reaction



Click here to watch a video



| Size | Knee Circumference | Thigh Circumference | SKU |
|---------|--------------------|---------------------|------------|
| Size XS | 30 - 34 cm | 40 - 44 cm | MGSK-XS-91 |
| Size S | 34 - 38 cm | 44 - 48 cm | MGSK-SM-91 |
| Size M | 38 - 42 cm | 48 - 52 cm | MGSK-MD-91 |
| Size L | 42 - 46 cm | 52 - 56 cm | MGSK-LG-91 |
| Size XL | 46 - 50 cm | 56 - 60 cm | MGSK-XL-91 |
| Size 2X | 50 - 54 cm | 60 - 64 cm | MGSK-2X-91 |

GENUACTIVE PAD

Patellar Knee Brace



Contains natural rubber latex that may cause allergic reaction



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




Compression | Stabilization | Massaging effect | Proprioception

FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap
- Anatomically shaped 3D knit for comfort and fit

COMPOSITION

52% Polyamide, 26% Thermoplastic Urethane, 14% Elastodiene, 5% Polyester, 2% Elastane, 1% Silicone



Measure knee circumference below the kneecap with 30 degrees of inclination, and thigh circumference 15 cm above the knee




Click here to watch a video

| Size | Knee Circumference | Thigh Circumference | SKU |
|---------|--------------------|---------------------|------------|
| Size XS | 30 - 34 cm | 40 - 44 cm | MGAK-XS-91 |
| Size S | 34 - 38 cm | 44 - 48 cm | MGAK-SM-91 |
| Size M | 38 - 42 cm | 48 - 52 cm | MGAK-MD-91 |
| Size L | 42 - 46 cm | 52 - 56 cm | MGAK-LG-91 |
| Size XL | 46 - 50 cm | 56 - 60 cm | MGAK-XL-91 |
| Size 2X | 50 - 54 cm | 60 - 64 cm | MGAK-2X-91 |


3D Knit
Provides breathability and comfort for both everyday and physical activities.

Lateral Supports
Provides high stability

Silicone Insert
The anatomical shape helps the kneecap to remain stable



Bend Zone
Knitted fabric is less compressive for more comfort and flexibility



Ankle

| | MALLEO CARE | MALLEO ACTIVE PAD | MALLEO ACTIVE BAND | MALLEO SUPPORT | MALLEO INTENSE |
|--|--------------------|----------------------|----------------------|----------------------|----------------------------|
| INDICATIONS | | | | | |
| Arthritis | ○ | ○ | ○ | ○ | ○ |
| Joint fluid and swelling | ○ | ○ | ○ | | |
| Chronic irritation/ Tendomyopathy | ○ | ○ | ○ | ○ | ○ |
| Post trauma, Post surgery, Soft tissue lesion | ○ | | | ○ | ○ |
| Hyperlaxity | | | | ○ | ○ |
| Chronic instability | | | | ○ | ○ |
| Post trauma, Post surgery, Ligament injury | | | | ○ | ○ |
| Moderate sprains | | | | ○ | ○ |
| SUPPORT LEVEL | Support level mild | Support level medium | Support level medium | Support level strong | Support level extra strong |

MALLEOCARE

Ankle Support



Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy, post trauma, post surgery, soft tissue lesion

ACTION MODE



Compression



Stabilization



Warming effect



Proprioception

FEATURES

- Mild support
- Adjustable level of stability
- Can be worn on either ankle

COMPOSITION

70% Polyamide, 12% Elastane, 9% Elastodiene, 8% Polyester, 1% Polypropylene



Measure ankle circumference above ankle bone



| Size | Circumference | SKU |
|---------------|---------------|------------|
| Size S | 18 - 22 cm | MMCA-SM-66 |
| Size M | 22 - 25 cm | MMCA-MD-66 |
| Size L | 25 - 28 cm | MMCA-LG-66 |

MALLEOACTIVE PAD

Ankle Brace



Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy

ACTION MODE



Compression



Stabilization



Massaging effect



Proprioception

FEATURES

- 2 lateral pads stabilize and support the ankle
- Gentle massaging effect
- Breathable 3D knit for optimal fit and comfort
- Left and right foot option
- Includes a foot slip for easier donning

COMPOSITION

45% Polyamide, 40% Thermoplastic Urethane, 10% Elastodiene, 3% Elastane, 2% Polyester



Measure ankle circumference above ankle bone



[Click here to watch a video](#)



| Size | Circumference | SKU Left Ankle | SKU Right Ankle |
|----------------|---------------|----------------|-----------------|
| Size S | 16 - 21 cm | MMAA-SM-L-91 | MMAA-SM-R-91 |
| Size M | 21 - 26 cm | MMAA-MD-L-91 | MMAA-MD-R-91 |
| Size L | 26 - 31 cm | MMAA-LG-L-91 | MMAA-LG-R-91 |
| Size XL | 31 - 37 cm | MMAA-XL-L-91 | MMAA-XL-R-91 |

MALLEOACTIVE BAND

Ankle Brace



Support level
medium



Contains natural rubber latex that may cause allergic reaction



[Click here to watch a video](#)



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy

ACTION MODE



Compression



Stabilization



Massaging effect



Proprioception

FEATURES

- The elastic band allows for adjustable support
- The anatomically shaped breathable knit fabric is compressive
- Can be worn in a shoe
- Can be worn on either ankle

COMPOSITION

62% Polyamide, 18% Elastodiene, 12% Viscose, 5% Polyester, 3% Elastane



Measure ankle circumference above ankle bone

| Size | Circumference | SKU |
|----------------|---------------|------------|
| Size XS | 18 - 20 cm | MMBA-XS-91 |
| Size S | 20 - 23 cm | MMBA-SM-91 |
| Size M | 23 - 25 cm | MMBA-MD-91 |
| Size L | 25 - 27 cm | MMBA-LG-91 |
| Size XL | 27 - 30 cm | MMBA-XL-91 |

MALLEOSUPPORT

Ankle Brace



Support level
strong

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate sprains

ACTION MODE



Movement control



Stabilization



Proprioception

FEATURES

- Stabilizing ankle brace with lateral rigid rods
- Two hook-and-loop straps for added support
- The lace closure system allows for an adjustable fit
- Can be worn on either ankle

COMPOSITION

91% Polyamide, 5% Elastane, 2% Polyester, 2% Cotton



Measure ankle circumference above ankle bone



| Size | Circumference | SKU |
|---------------|---------------|------------|
| Size S | 19 - 22 cm | MMSA-SM-91 |
| Size M | 22 - 26 cm | MMSA-MD-91 |
| Size L | 26 - 30 cm | MMSA-LG-91 |

MALLEOINTENSE

Ankle Stabilizer



Contains natural rubber latex that may cause allergic reaction




Click here to watch a video



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate to severe sprains

ACTION MODE




FEATURES

- Immobilizing ankle stabilizer based on two rigid shells
- Adaptable fit with two gel cushions to provide comfort and cooling effect
- Two Velcro straps as well as heel support strap provide ability to adjust tightness and support
- Can be worn on either ankle

COMPOSITION

50% Gel, 22% Polyvinyl chloride, 20% Polypropylene, 7% Polyamide, 1% Neoprene



| Size | Height | SKU |
|----------|----------|---------|
| One Size | ≥ 150 cm | MMIA-91 |

HOW TO ORDER MOBILIS PRODUCTS

| M | E | A | E | - | MD | - | - | - | 91 |
|---------|--|--|---|---|------|---|--|---|--|
| Mobilis | Letter matches first part of the product name: E = Epi M = Manu MP = Manu Plus R = Rhizo L = Lumbo G = Genu M = Malleo | Letter matches line of product: C = Care W = Wrap A = Active S = Support or Stab B = Band I = Intense SE = Support Evo | Letter matches body part where the product is used: E = Elbow W = Wrist T = Thumb B = Back K = Knee A = Ankle | | Size | | When your product has a left or right side, L or R is indicated. When your product has a short or tall height, S or T is indicated. | | Available colours: 66 = Beige 91 = Graphite 99 = Black |

Add **L** to indicate the Left side or **R** to indicate the Right side after the Size in the ordering code. **L** or **R** is not specified in products that fit both left and right. **Example:** MMAW-MD-L-91.

For back braces, add **S** to indicate the Short height or **T** to indicate the Tall height after the Size in the ordering code.

You can find products without a size. **Example:** MEBE-91

You can find products with combined sizes. **Example:** MRST-SM/MD-91